Quick access to neutralizing eyewash is crucial in case of chemical accidents

Plum Eyewash pH Neutral 4.9% Phosphate is a buffered solution that quickly neutralizes acids and alkali

PRODUCT FACTS & BENEFITS

- Quick neutralising effect - minimizes the impact of accidents with acids and alkali
- Easy to handle - easily applied by the injured person or with assistance
- Flexible solutions - bottles for mobile use and stations for wallmounting
- Wide application - from industrial workplaces to laboratories and first aid kits
- High functional reliability - bottles comply with to EN 15154-4
- CE marked according to the Guidelines for Classification of MDD 93/42/EEC
- Shelf life of 3.5 years for sealed bottles

IMPORTANT

Reservation is made for hydrofluoric acid, which is toxic and highly corrosive to eyes, skin, and bones. Plum Eyewash pH Neutral 4.9% Phosphate is not intended as first aid for eye accidents involving hydrofluoric acid, because it does not limit the corrosive effect of HF.

info@plum.eu
+45 6916 9600

Optimize safety with a quick, effective and accessible eyewash solution


Using the Plum Eyewash pH Neutral 4.9% Phosphate buffer quickly neutralises the pH value of the chemical to a harmless level.

A quick neutralising effect is key in minimising the impact of an accident with acids and alkali.

Flushing the eye with normal water or Plum Eyewash 0.9% Sodium Chloride solution only creates a diluting and rinsing effect (see figure 1).

After using Plum Eyewash pH Neutral 4.9% Phosphate continue rinsing with Plum Eyewash 0.9% Sodium Chloride until you reach a doctor.

Figure 1: Neutralizing effect of Plum Eyewash pH Neutral 4.9% Phosphate and 0.9% sodium chloride

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ERGONOMIC EYECUPS

The bottles are equipped with a selection of ergonomic eyecups to ensure that the liquid can be applied correctly and in an even flow. Squeeze the bottle gently to generate a soft and pleasant flow.

EYECUP

Flexible, ergonomic eyecup that keeps the eye lid open during use

DUO EYECUP

Flexible, ergonomic eyecup that allows rinsing of both eyes at the same time

AVAILABLE SIZES

Art. no.	4752	4755 - single carton	4801
Size	200 ml	200 ml	500 ml
Packaging	10 x 200 ml	10 x 200 ml	6 x 500 ml

info@plum.eu
+45 6916 9600