

# Safety Data Sheet

## Aloe Vera Burn Gel

Issued: 18/02/2019  
Version: 1.0.0

### SECTION 1: Identification of the substance/preparation and of the company/undertaking

#### 1.1. Product identifier

**Trade name:** Aloe Vera Burn Gel

**Other Information:** The product is a medical device covered by Council Directive 93/42/EEC on medical devices. Consequently, the product does not require a safety data sheet and therefore this data sheet does not fully meet the statutory requirements for a safety data sheet.

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

**Recommended uses:** This safety data sheet contains general information concerning the chemical, but specific instructions and guidelines are indicated on the product label and in the instructions for use.

#### 1.3. Details of the supplier of the safety data sheet

##### Manufacturer

**Company:** Aurena Laboratories AB

**Address:** Fjärrviksvägen 22

**Zip code:** 653 50

**City:** Karlstad

**Country:** SWEDEN

##### Distributor

**Company:** Plum A/S

**Address:** Frederiks Plums Vej 2

**Zip code:** Dk-5610

**City:** Assens

**Country:** DENMARK

#### 1.4. Emergency Telephone Number

Members of the public: 111 (NHS 111 (Scotland: NHS 24)) The emergency telephone is open 24 hours.

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

**CLP-classification:** Aerosol 3;H229

**Most serious harmful effects:** Pressurised container: May burst if heated.

#### 2.2. Label elements

**Signal word:** Warning

##### H-phrases

H229 Pressurised container: May burst if heated.

##### P-phrases

P102 Keep out of reach of children.

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

# Safety Data Sheet

## Aloe Vera Burn Gel

Issued: 18/02/2019  
Version: 1.0.0

P251	Do not pierce or burn, even after use.
P410	Protect from sunlight.
P412	Do not expose to temperatures exceeding 50°C/122°F.

### 2.3. Other hazards

The product does not contain any PBT or vPvB substances.

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures

Contains no substances subject to reporting requirements

**Ingredient comments:** Composition: Purified water, Carbomer, Preservatives, Aloe Vera, Sodium Hydroxide.

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

<b>Inhalation:</b>	Seek fresh air.
<b>Ingestion:</b>	Wash out mouth thoroughly and drink 1-2 glasses of water in small sips. Seek medical advice in case of persistent discomfort.
<b>Skin contact:</b>	In case of irritation: Wash the skin thoroughly with water and continue washing for a long time. Seek medical advice in case of persistent discomfort.
<b>Eye contact:</b>	Flush with water (preferably using eye wash equipment) until irritation subsides. Seek medical advice if symptoms persist.
<b>General:</b>	When obtaining medical advice, show the safety data sheet or label.

### 4.2. Most important symptoms and effects, both acute and delayed

May cause eye irritation.

### 4.3. Indication of any immediate medical attention and special treatment needed

No special immediate treatment required. Treat symptoms.

## SECTION 5: Fire-fighting measures

### 5.1. Extinguishing media

**Suitable extinguishing media:** Extinguish with powder, foam, carbon dioxide or water mist. Use water or water mist to cool non-ignited stock.

**Unsuitable extinguishing media:** Do not use water stream, as it may spread the fire.

### 5.2. Special hazards arising from the substance or mixture

CAUTION! Aerosol containers may explode.  
Heating will cause a rise in pressure in packaging with a risk of bursting. Can generate harmful flue gases containing carbon monoxide in the event of fire.

### 5.3. Advice for fire-fighters

If there is a risk of exposure to vapour and flue gases, a self-contained breathing apparatus must be worn.  
Move containers from danger area if it can be done without risk. Avoid inhalation of vapour and flue gases - seek fresh air.

## SECTION 6: Accidental release measures

# Safety Data Sheet

## Aloe Vera Burn Gel

Issued: 18/02/2019  
Version: 1.0.0

### 6.1. Personal precautions, protective equipment and emergency procedures

**For non-emergency personnel:** Keep unnecessary personnel away. Smoking and naked flames prohibited.

**For emergency responders:** In addition to the above: Normal protective clothing equivalent to EN 469 is recommended.

### 6.2. Environmental precautions

Avoid unnecessary release to the environment.

### 6.3. Methods and material for containment and cleaning up

Contain and absorb spill with sand or other absorbent material and transfer to suitable waste containers. Wipe up minor spills with a cloth.

### 6.4. Reference to other sections

See section 8 for type of protective equipment. See section 13 for instructions on disposal.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

Smoking and naked flames prohibited. Running water and eye wash equipment should be available. Wash hands before breaks, before using restroom facilities, and at the end of work.

### 7.2. Conditions for safe storage, including any incompatibilities

Store safely and keep out of reach of children. Pressurized container: Protect from sunlight and do not expose to temperatures exceeding 50°C. Storage temperature: 15-30 °C.

### 7.3. Specific end use(s)

None.

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

**Occupational exposure limit:** Contains no substances subject to reporting requirements

**Legal basis:** EH40/2005 Workplace exposure limits. Last amended December 2011.

**Other Information:** The product requires no special exposure control.

### 8.2. Exposure controls

**Appropriate engineering controls:** Wear the personal protective equipment specified below.

**Personal protective equipment, eye/face protection:** Not required.

**Personal protective equipment, skin protection:** Not required.

**Personal protective equipment, hand protection:** Not required.

**Personal protective equipment, respiratory protection:** Not required.

**Environmental exposure controls:** Ensure compliance with local regulations for emissions.

# Safety Data Sheet

## Aloe Vera Burn Gel

Issued: 18/02/2019  
Version: 1.0.0

### SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Parameter	Value/unit
State	Aerosol Gel
Colour	Clear Colourless
Odour	Neutral
Solubility	Soluble in the following: Water.
Explosive properties	Non-explosive
Oxidising properties	Test data are not available.

Parameter	Value/unit	Remarks
pH (solution for use)	No data	
pH (concentrate)	6.0 - 7.5	
Melting point	No data	
Freezing point	No data	
Initial boiling point and boiling range	~ 100 °C	
Flash Point	> 65 °C	
Evaporation rate	No data	
Flammability (solid, gas)	No data	
Flammability limits	No data	
Explosion limits	No data	
Vapour pressure	No data	
Vapour density	No data	
Relative density	No data	
Partition coefficient n-octanol/water	No data	
Auto-ignition temperature	No data	
Decomposition temperature	No data	
Viscosity	25000 - 60000 cps	
Odour threshold	No data	

#### 9.2 Other information

**Other Information:** None.

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

Not reactive.

#### 10.2. Chemical stability

The product is stable when used in accordance with the supplier's directions.

#### 10.3. Possibility of hazardous reactions

No risk of hazardous reactions.

#### 10.4. Conditions to avoid

Do not expose to heat (e.g. sunlight).

#### 10.5. Incompatible materials

## Safety Data Sheet

### Aloe Vera Burn Gel

Issued: 18/02/2019  
Version: 1.0.0

None known.

#### 10.6. Hazardous decomposition products

None known.

### SECTION 11: Toxicological information

#### 11.1. Information on toxicological effects

<b>Acute toxicity - oral:</b>	Ingestion may cause discomfort. The product does not have to be classified.
<b>Acute toxicity - dermal:</b>	The product does not have to be classified.
<b>Acute toxicity - inhalation:</b>	The product does not have to be classified.
<b>Skin corrosion/irritation:</b>	The product does not have to be classified.
<b>Serious eye damage/eye irritation:</b>	May cause eye irritation. The product does not have to be classified.
<b>Respiratory sensitisation or skin sensitisation:</b>	The product does not have to be classified.
<b>Germ cell mutagenicity:</b>	The product does not have to be classified.
<b>Carcinogenic properties:</b>	The product does not have to be classified.
<b>Reproductive toxicity:</b>	The product does not have to be classified.
<b>Single STOT exposure:</b>	The product does not have to be classified. Test data are not available.
<b>Repeated STOT exposure:</b>	The product does not have to be classified. Test data are not available.
<b>Aspiration hazard:</b>	The product does not have to be classified.

### SECTION 12: Ecological information

#### 12.1. Toxicity

The product does not have to be classified.

#### 12.2. Persistence and degradability

Test data are not available.

#### 12.3. Bioaccumulative potential

No bioaccumulation expected.

#### 12.4. Mobility in soil

Test data are not available.

#### 12.5. Results of PBT and vPvB assessment

The product does not contain any PBT or vPvB substances.

#### 12.6. Other adverse effects

None known.

# Safety Data Sheet

## Aloe Vera Burn Gel

Issued: 18/02/2019  
Version: 1.0.0

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

Do not dispose of aerosol sprays in refuse collection, even when empty. The sprays must be sent to the municipal chemical waste collection facility with the specifications set out below.

**Category of waste:** 16 05 05 gases in pressure containers other than those mentioned in 16 05 04

### SECTION 14: Transport information

#### Land transport (ADR/RID)

<b>14.1. UN-No.:</b> 1950 <b>14.2. UN proper shipping name:</b> AEROSOLS  <b>14.3. Transport hazard class(es):</b> 2.2 <b>Hazard label(s):</b> 2.2 <b>Hazard identification number:</b>	<b>14.4. Packing group:</b> <b>14.5. Environmental hazards:</b> The product should not be labelled as an environmental hazard (symbol: fish and tree).  <b>Tunnel restriction code:</b> E
--	--

#### Inland water ways transport (ADN)

<b>14.1. UN-No.:</b> 1950 <b>14.2. UN proper shipping name:</b> AEROSOLS  <b>14.3. Transport hazard class(es):</b> 2.2 <b>Hazard label(s):</b> 2.2 <b>Transport in tank vessels:</b>	<b>14.4. Packing group:</b> <b>14.5. Environmental hazards:</b> The product should not be labelled as an environmental hazard (symbol: fish and tree).
---	---

#### Sea transport (IMDG)

<b>14.1. UN-No.:</b> 1950 <b>14.2. UN proper shipping name:</b> AEROSOLS <b>14.3. Transport hazard class(es):</b> 2.2 <b>Hazard label(s):</b> 2.2 <b>EmS:</b> F-D, S-U	<b>14.4. Packing group:</b> <b>14.5. Environmental hazards:</b> The product is not a Marine Pollutant (MP). <b>Environmental Hazardous Substance Name(s):</b> <b>IMDG Code segregation group:</b> - None -
--	---

#### Air transport (ICAO-TI / IATA-DGR)

<b>14.1. UN-No.:</b> 1950 <b>14.2. UN proper shipping name:</b> AEROSOLS, NON-FLAMMABLE <b>14.3. Transport hazard class(es):</b> 2.2 <b>Hazard label(s):</b> 2.2	<b>14.4. Packing group:</b> <b>14.5. Environmental hazards:</b>
---	--

#### 14.6. Special precautions for user

None.

#### 14.7. Transport in bulk according to Annex II of MARPOL and the IBC Code

Not applicable.

# Safety Data Sheet

## Aloe Vera Burn Gel

Issued: 18/02/2019  
Version: 1.0.0

**Other Information:** May be transported in limited quantities if placed in outer packaging according to ADR 3.4, when max. 1 liter/inner packaging and max. 30 kg total gross mass. Shrink- or stretch wrapped trays may be used and shall not exceed 20 kg total gross mass/tray.  
May be transported in limited quantities if placed in outer packaging according to IMDG 3.4, when max. 1 liter/inner packaging and max. 30 kg total gross mass. Shrink- or stretch wrapped trays may be used and shall not exceed 20 kg total gross mass/tray.

### SECTION 15: Regulatory information

#### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

**Special Provisions:** The product is covered by Council Directive 93/42/EEC on medical devices. Consequently, the product is not covered by Regulation (EC) no 1272/2008 of the European Parliament and the Council on classification, labelling and packaging of substances and mixtures. cf. Article 1 (5d). Regulation 1907/2006 EC of the European Parliament and Council on registration, evaluation, authorisation and restriction of chemicals (REACH) states in Article 2 (6c) that section IV regarding requirements for safety data sheets does not apply to medical equipment. Consequently, there is no requirement stating that the product mentioned above should be provided with a safety data sheet to meet the requirements of the above regulation.

European Agreement concerning the International Carriage of Dangerous Goods by Road (UN ADR 2019).  
COUNCIL DIRECTIVE of 20 May 1975 on the approximation of the laws of the Member States relating to aerosol dispensers (75/324/EEC)

#### 15.2. Chemical Safety Assessment

**Other Information:** Chemical safety assessment has not been performed.

### SECTION 16: Other information

#### Version history and indication of changes

Version	Revision date	Responsible	Changes
1.0.0	18/02/2019	CCS Hygien A/S (SKP)	New

**Abbreviations:**  
DNEL: Derived No Effect Level  
PNEC: Predicted No Effect Concentration  
PBT: Persistent, Bioaccumulative and Toxic  
vPvB: Very Persistent and Very Bioaccumulative  
STOT: Specific Target Organ Toxicity

**Other Information:** Expiry date for unopened, sealed bottles: At least three years from the manufacturing date. Additional information can be obtained by contacting info@plum.dk or by telephone +45 6471 2112.

**Date:** 18/02/2019

#### List of relevant H-statements

H229 Pressurised container: May burst if heated.

**Product type:** MDD Class II a according to MDD 93/42 EEC.

#### SDS is prepared by

**Company:** Plum A/S

**Address:** Frederik Plums Vej 2

## Safety Data Sheet

### Aloe Vera Burn Gel

Issued: 18/02/2019  
Version: 1.0.0

**Zip code:** 5610  
**City:** Assens  
**Country:** DENMARK  
**E-mail:** info@plum.dk  
**Phone:** +45 64712112  
**Fax:** +45 64712125  
**Homepage:** www.plum.dk  
**Document language:** GB