


QuickFix® plaster refills


Plum has launched 3 new refills for our popular QuickFix® plaster dispenser

For our range of plaster solutions we introduce 3 new refills: Aloe Vera, Alu and Blood Stopper. Each plaster is packed individually and ready for use when pulled from either the original QuickFix® dispenser or the new QuickFix® UNO dispenser. Each plaster can be affixed with only one hand resulting in an easy and quick return back to work.



INTENDED USE: Small cuts and abrasions

DIRECTIONS: Apply on well cleaned, dry skin



QUICKFIX® ALOE VERA REFILL

Article no: 5514
Packaging: 6 x 45 pcs.
Dimensions: L: 72 x W: 25 mm

The aloe vera contained in the pad has well known soothing, moisturizing and calming properties. It creates the best conditions on the wound environment helping the natural healing process. Prevents the formation of scabs.

PRODUCT CLAIMS:

- Soothing and moisturizing
- Helps the healing process
- Air-permeable
- Skin friendly
- Suitable for sensitive skin




QUICKFIX® ALU REFILL

Article no: 5515
Packaging: 6 x 45 pcs.
Dimensions: L: 72 x W: 25 mm

Plaster with micronized aluminium pad. The anti-adhering layer in aluminium and PE has a higher anti-adhering performance compared to a standard PE layer. The plaster is air-permeable and comfortable.

PRODUCT CLAIMS:

- Improved anti-adhering effect
- Air-permeable
- Skin friendly
- Suitable for sensitive skin



QUICKFIX® BLOOD STOPPER REFILL

Article no: 5516
Packaging: 6 x 45 pcs.
Dimensions: L: 72 x W: 25 mm

Plasters with haemostatic effect. The presence of alginate fibres between the nonwoven fibres of the wound pad creates an artificial clot (gel) inside the pad, which blocks the bleeding in 1-2 minutes and maintains optimal moist conditions.

PRODUCT CLAIMS:

- Haemostatic plasters
- Alginate fibers
- Improved protection effect
- Helps the healing process

DISPENSER SOLUTIONS:

The new QuickFix® refills are available as dispenser solutions with 2 refills in the QuickFix® dispenser.

They also fit in the QuickFix® UNO dispenser with either a white or transparent cover.



5505 QuickFix® ALU
 5506 QuickFix® ALOE VERA
 5510 QuickFix® BLOOD STOPPER



5534 QuickFix® UNO, White
 5533 QuickFix® UNO, Transparent