Technical Data Sheet Plum Eyewash 0.9% Sodium Chloride

1. Identification of the product and the company

Trade Name:	Plum Eyewash 0.9% Sodium Chloride and Plum Eyewash 0.9% Sodium Chloride QuickRinse		
Intended Use:	The products are intended to be used as first aid products for eye irrigation.		
Expiry date:	Printed on the bottle label/ampoule.		
1.2. Details of the supplier o	f the data sheet		
Manufacturer			
Company:	Plum Safety ApS		
Address:	Mandelalleen 1		
Zip code:	5610		
City:	Assens		
Country	DENMARK		
E-mail	info@plum.eu		
Phone	+45 6916 9600		
Fax	+45 6916 9630		
1.3 Legislative information General information:	As the products are class 1s sterile medic		
	 93/42/EEC of 14 June 1993 concerning medi Body, DNV Product Assurance AS. No devices. The products shall not be classified as hazard labeling rules for substance and mixtures. Plum Safety ApS is certified according to ISG management systems – Requirements for regulatory put 	Safety Data Sheet is required for medica lous according to the classification and <i>O</i> 13485:2016 Medical devices-Quality	
1.4 Instruction for use:	Mono eye cup (200, 500, 1000 ml) 1 2 3 1 360° 2 3	Duo eye cup (500, 1000 ml) $ \begin{array}{c} 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\$	
	Ampoule (20 ml)		
	- Do not reuse.	\otimes	
	- Do not use if damaged.		
	- Recommended temperature when used	16°C	

Technical Data Sheet Plum Eyewash 0.9% Sodium Chloride

2. General information about the product

Material requirement	its:	The materials in these products have been judged biocompatible according the		
		requirements of ISO 10993-1:2009 Biological evaluation of medical devices - Part 1:		
		Evaluation and testing within a risk management process.		
Bottle, label and eye cup dimensions.				
	20 ml ampoule	200 ml bottle	500 ml bottle	1000 ml bottle

, ,				
	20 ml ampoule	200 ml bottle	500 ml bottle	1000 ml bottle
Bottle material	Polyethylene	Polyethylene without	Polyethylene without	Polyethylene without
	without additives	additives	additives	additives
The bottle/ampoule	20.0-20.5	201-207	505-512	1005-1020
contains [ml]				
Bottom Ø [mm]	21	58.5	78	90
Height [mm]	104	129.6	183	229

3. Side effects and malfunction		
Side effects and malfunction:	The products are designed as a first aid product for eye irrigation. If side effects or malfunctions are experienced, please report these to Plum Safety ApS, info@plum.eu or	
	phone +45 6916 9600, so that necessary remedial measures can be initiated.	
	Serious adverse effects should always be reported to a relevant competent authority.	

4. Storage		
Storage:	5°C	
Shelf life:	42 months from manufacturing date.	
Installation:	Please refer to the pictogramm on the label:	
Mounting stations and boxes:	Please visit https://plum.eu/products/firstaid/eye-wash	

5. Other information

Version history and indication of changes			
Version	Revision date	Responsible	Changes
1.3	25-05-2021	CL	1.2, 1.3, 3, 5
1.2	30-09-2020	CL	5
1.1	03-06-2020	CL / VWB	5
1.0	14-04-2020	CL /VWB	

Other information Additional information can be obtained by contacting <u>info@plum.eu</u> or by telephone +45 6916 9600

TDS prepared by:	
Company:	Plum Safety ApS
Address:	Mandelalleen 1
Zip code:	5610
City	Assens
Country:	DENMARK
E-mail:	info@plum.eu
Phone:	+45 6916 9600
Fax:	+45 6916 9630
Document Language:	English