# **Technical Data Sheet**

## **OuickFix**

(elastic, detectable, water resistant & alu)

## 1. Identification of the product and the company

#### 1.1. Product identifier and intended use.

Trade Name: QuickFix and QuickFix Long

**Intended Use:** Plum QuickFix products are intended for use as a mechanical barrier for absorption

of exudates from minor wounds.

**Expiry date:** Printed on the back of the plaster.

#### 1.2. Details of the supplier of the data sheet

Manufacturer

Company: Plum Safety ApS
Address: Mandelalleen 1

 Zip code:
 5610

 City:
 Assens

 Country
 DENMARK

 E-mail
 info@plum.eu

 Phone
 +45 6916 9600

 Fax
 +45 6916 9630

#### 1.3 Legislative information

General information: As the products are class 1 medical devices according to Council Directive

93/42/EEC of 14 June 1993 concerning medical devices. No notified body is

involved. No Safety Data Sheet is required for medical devices.

Plum Safety ApS is certified according to ISO 13485:2016 Medical devices - Quality

management systems - Requirements for regulatory purposes.

## 2. General information about the product

Material requirements:

The materials in these products have been judged biocompatible according to the requirements of ISO 10993-1:2009 Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process.

#### Physical properties and composition:

	elastic	detectable	water resistant	alu
Material	Pink cross direction elastic fabric (Cotton & polyamide)	Blue cross direction elastic fabric (Cotton & polyamide)	Pink perforated PE (Polyethylene)	Pink cross direction fabric (Cotton & polyamide)
Adhesive pad	Hypoallergenic acrylic polymer	Hypoallergenic acrylic polymer	Hypoallergenic acrylic polymer	Hypoallergenic acrylic polymer
Absorbing pad	Viscose and non-woven PP	Viscose and PP/PE non-woven	Viscose and non-woven PP	Viscose and aluminium coated PE non-woven
Primary Packing (Pouch)	Cold sealed paper	Cold sealed paper	Cold sealed paper	Cold sealed paper
Dimensions	72 x 25 mm 120 x 20 mm	72 x 25 mm 120 x 20 mm	72 x 25 mm	72 x 25 mm
Special feature	Elastic	Elastic Detectable aluminium foil facing the backing	Water resistant	Elastic Anti-adhering effect

# **Technical Data Sheet**

## **OuickFix**

(elastic, detectable, water resistant & alu)

## 3. Side effects and malfunction

Side effects and malfunction: The products are designed for skin contact and first aid use involving minor wounds.

If side effects or malfunctions are experienced, please report these to Plum Safety ApS, <a href="mailto:info@plum.eu">info@plum.eu</a> or phone +45 6471 2112, so that necessary remedial measures can be

initiated.

Serious adverse effects should always be reported to a relevant and competent

authority.

## 4. Storage

Storage:





Keep dry

Keep away from sunlight

**Shelf life:** 5 years from production date.

## 5. Other information

## Version history and indication of changes

Version	Revision date	Responsible	Changes
1.3	25-05-2021	CL	1.2, 1.3, 3, 5
1.2	30-09-2020	CL	5
1.1	12-06-2020	CL	5
1.0	14-04-2020	CL	

Other information Additional information can be obtained by contacting info@plum.eu or by telephone +45 6471 2112.

TDS prepared by:

**Company:** Plum Safety ApS **Address:** Mandelalleen 1

 Zip code:
 5610

 City
 Assens

 Country:
 DENMARK

 E-mail:
 info@plum.eu

 Phone:
 +45 6916 9600

 Fax:
 +45 6916 9630

**Document Language:** English