



SAFETY AT WORK

CUSTOMIZED
FIRST AID
SOLUTIONS

plum
Safety

Welcome to Plum Safety

We are international experts and offer high quality products and expert knowledge within first aid solutions and hygiene for the workplace.

We want you to be safe at work!

WORKPLACE HYGIENE & FIRST AID

Accidents happen every day. Minimize the extent and impact by making sure you have the right first aid solutions available where accidents can happen.

A first aid solution tailored to your workplace will help you prepare for accidents and quickly treat eye injuries, cuts and minor wounds in a safe and hygienic way.

A wide range of eyewash, plasters, decontamination solutions, burn gel, and wound cleansing wipes, helps to quickly treat you or your colleague and minimize the damage of an accident.

It is essential that the mix of safety products is fitted to the requirements of your workplace.

GLOBAL FOOTPRINT

We have a global footprint and are represented in most of the world. Reach out if you need us to help you find your safety supplier!



You are welcome to contact us
info@plum.eu
+45 69 16 96 00

Table of content


Protect your people	4
Plum Safety Reminder App	5
How to choose the correct eyewash	6
Eyewash	8
pH Neutral	10
Wound & Eyewash Spray	12
Open eyewash stations	14
Plum iBox	17
Eyewash boxes	20
Decontamination solution - Plum BodyNeutrAll	22
QuickFix plasters & dispensers	24
QuickSafe®.....	30
Plum Pull1Aid®	34
Wound care	36
Fire Safety	38
Plum Safety in brief	39



For the online version: click on the link and jump directly to the selected page.

Protect your people Equip your facility right

Quick and easy access to first aid will limit the impact of accidents and keep your workplace safe



Accidents happen every day. Minimize the extent and impact by making sure you have the right first aid solutions available where accidents can happen. A first aid solution tailored to your workplace will help you prepare for accidents and quickly treat eye injuries, cuts and minor wounds in a safe and hygienic way.

QUICK AND EFFECTIVE HELP TAILORED TO YOUR BUSINESS

A wide range of eyewash, plasters, burn gel and wound cleansing wipes, helps to quickly treat you or your colleague and minimize the damage of an accident. It is essential that the mix of safety products is fitted to the requirements of your workplace.

We will help you find the right solution, whether you work in an office or handle acids or alkali in a high-risk environment like charging batteries, industrial cleaning, laboratory work or petrochemical industry.

PRACTICAL SOLUTIONS FOR BUSY WORKING DAYS

Your hands, your eyes and your head are your most important working tools. Even a small incident can affect your work, and larger ones may cause an even greater impact if you are not prepared. By taking the correct precautions, you enable yourself to have the right first aid within reach, whereby many smaller everyday accidents can be minimized significantly.



Plum Eyewash is easy to use if an accident happens



Turn the eyecup in the direction of the arrow 360° until the seal breaks



Lean back, squeeze the bottle gently and rinse



Or lean forward and rinse by squeezing the bottle gently

Free account. Unlimited products. Auto notifications. Log inspections.

SIMPLE & EFFICIENT TOOL FOR MANAGEMENT OF EXPIRY DATES


Help your customer and position yourself as dealer in the app Plum Safety Reminder - a free tool for efficient management of expiry dates

Developed to make life easier for employees responsible for safety products by offering a complete product overview and automatic email alerts.

FEATURES

- ✓ Free account
- ✓ Simple management of expiration dates
- ✓ Automatic email alerts
- ✓ Manage multiple locations
- ✓ ... and multiple dealers
- ✓ Scan products to load information
- ✓ Log function for inspections

+45 6916 9600
info@plum.eu



**SCAN &
DOWNLOAD
NOW**



plum
Safety

How to choose the correct eyewash

When it comes to eye accidents, the harmful substance varies and you need to choose the right eyewash accordingly

The cause of an accident determines which eyewash is required. When an accident occurs, it is essential to determine whether it is caused by either a foreign object or by acids or alkali - and if possible, which pH value the substance in question holds.

Know your chemicals and ensure the best suited first aid.

An eyewash plan helps you choose the right eyewash by clear guidelines. Place the eyewash plan next to the eyewash stations, boxes and QuickSafe®.



How to use Plum Eyewash



How to use Plum Eyewash DUO



Art.nr. 956480B

Dimensions: W: 22,8 x H: 29,0 cm

The pH scale ranges from 0 to 14 and measures how acidic or alkalic a substance is.

- A pH of 7 is neutral
- A pH less than 7 is acidic
- A pH greater than 7 is alkalic

When pH is either strongly acidic or alkalic, start by rinsing with Plum pH Neutral until the bottle is empty. Discard the bottle and rinse with Plum Eyewash until you reach a doctor or the emergency room.

When pH is neutral, rinse with Plum Eyewash until the bottle is empty and discard the bottle after use.

QR-code :
How to use Plum Eyewash :
(animation video) :

QR-code :
How to use Plum Eyewash DUO :
(animation video) :

Mount the eyewash plan on the wall :
next to the eyewash station. Mark off :
145 cm from the floor and place the :
eyewash plan according to the mark.



Rinse both eyes simultaneously

You no longer have to choose which eye to save. Plum DUO eyecup is a unique solution and allows you to rinse both eyes at the same time and save them both

MANY ACCIDENTS AFFECT BOTH EYES

Eye accidents involving chemicals often affect both eyes. A DUO eyecup helps you prepare, because it allows rinsing of both eyes at the same time. Stop worrying about which eye to save, save them both at once instead.

EYECUP FOR OPTIMAL RINSING

Apart from the strength of the acid or alkali, the contact time is the most critical factor in accidents involving acids and alkali. It is key that the available first aid neutralize the harmful substance as quickly as possible. With the unique DUO eyecup, you can treat both eyes at the same time.

The eyecup is shaped to help you keep your eyes open during rinsing and ensure the correct flow of liquid and rinsing time.

Plum Eyewash and Plum pH Neutral are disposable products, which ensure that the rinsing liquid is sterile every time



Eyewash in case of foreign objects

Plum Eyewash eases discomfort and prevents following complications caused by foreign bodies in the eyes

You can easily get a foreign object in your eye. Dirt, dust, cleaning agents or paint from your hands all pose a risk for your eyes. The risk increases when working with machines and tools that can cause wood or metal splinters in your eyes.

Getting a foreign body in your eye can be extremely unpleasant, as the eye is highly sensitive. Without treatment, it can develop into an inflammation of the eye, cornea or iris. By rinsing immediately after the incident, you can prevent this and decrease the risk of following complications.

Plum Eyewash contains a sterile saline solution of 0.9% sodium chloride, which is similar to the eye's natural fluid. Squeeze the bottle gently to generate a soft and pleasant flow. The eyewash product range includes everything from a small handy bottle that you can carry with you, to wall-mounted stations that makes it quick and easy to get help in case of an emergency.



**200 ML / 500 ML / 1000 ML
PLUM EYEWASH**

Bottles with a 0.9% sterile sodium chloride solution. The bottles can be used individually or as refills in eyewash stations and boxes. They come with a dust cap. The 200 ml bottle is easy to store in the specially designed belt bag, in the first aid kit, in the toolbox, etc.

Content: 200 ml
Packaging: 10 x 200 ml
Article no: 4691

Content: 500 ml
Packaging: 12 x 500 ml
Article no: 4604

Content: 1000 ml
Packaging: 6 x 1000 ml
Article no: 4707



**5 X 20 ML
QUICKRINSE REFILL**

5 eyewash ampoules with 20 ml 0.9% sterile sodium chloride solution. For use individually or for refilling QuickSafe®, QuickFix&Rinse and QuickRinse dispensers.

Content: 5 x 20 ml
Packaging: 4 x 5 x 20 ml
Article no: 5160



**500 ML / 1000 ML
PLUM EYEWASH DUO**

Bottles with a 0.9% sterile sodium chloride solution. The ergonomic DUO eyecup makes it possible to rinse both eyes at the same time.

Content: 500 ml
Packaging: 6 x 500 ml
Article no: 4861



**200 ML / 500 ML
PLUM EYEWASH SINGLE CARTONS**

Single cartons with Plum Eyewash. The carton protects the bottle. Suitable for mobile workplaces or for refilling stations.

Content: 200 ml
Packaging: 10 x 200 ml
Article no: 4695



**3 X 500 ML
PLUM EYEWASH ON-THE-GO**

Practical and handy set with 3 bottles of Plum Eyewash. Easy to take with you in the car or on the way to a doctor or the emergency room.

Content: 3 x 500 ml
Packaging: 4 x 3 x 500 ml
Article no: 4620

Content: 1000 ml
Packaging: 6 x 1000 ml
Article no: 4800

Content: 500 ml
Packaging: 12 x 500 ml
Article no: 4605



Eyewash in case of chemical accidents

Plum pH Neutral
neutralizes acids and
alkali and quickly flush
them from the eye
during rinsing

IMPORTANT



HF

Reservation is made for hydrofluoric acid, which is toxic and highly corrosive to skin and bones.

Plum pH Neutral and Plum Eyewash may neutralize and rinse, but not limit the corrosive effect.

They are not intended as first aid for eye accidents involving hydrofluoric acid.

RECOMMENDED USE OF PLUM PH NEUTRAL

Plum pH Neutral is an eyewash solution for accidents involving acids and alkali.

Rinse with Plum pH Neutral to reduce the pH value of the harmful substance to a harmless level as quickly as possible.

Continue rinsing with Plum Eyewash 0.9% sodium chloride until you reach a doctor or emergency room.

DANGEROUS CHEMICALS ARE EVERYWHERE

In many jobs, businesses and industries, there is a risk of being exposed to serious eye injuries caused by chemicals e.g. in chemical factories, the recycling industry, on construction sites and in the mining industry. You need a neutralizing eyewash to minimize the impact of the accident.

Acids and alkali can break down the structure of skin and mucosa membranes. The severity of corrosion injuries depends primarily on:

- The concentration of the chemical
Strong acids and alkali are extremely corrosive
- The temperature
A high temperature increases the corrosive properties of the chemical
- The duration of contact
The longer the duration, the deeper the harmful substance will penetrate

A FEW SECONDS CAN MAKE THE DIFFERENCE

In case of eye injuries involving acids and alkali, it is crucial to neutralize the chemical as quickly as possible to a pH level of approximately 7.4, corresponding to the normal pH level of the eye's fluid - before the acid or alkali has penetrated the cornea.

To optimize safety where there is a risk of acid and alkali injuries, the eyewash solution must be:


- Easily accessible
- Quick to use
- Effective

QUICK AND EFFECTIVE EYEWASH IS KEY

Plum pH Neutral is a 4.9% phosphate buffer that will quickly and effectively neutralize even strong acids and alkali. The liquid consists of phosphate salts that are present in the body. The hydrogen phosphate buffer is effective against acids, whereas the dihydrogen phosphate buffer reacts against alkali. The effect of Plum pH Neutral has been documented in scientific studies.^{1) 2)}

We have tested the neutralizing effect of Plum pH Neutral for both strong acids and alkali compared with the effect of saline eyewash with 0.9% sodium chloride. The test clearly shows that Plum pH Neutral will quickly neutralize both acids and alkali to a harmless pH level. Saline eyewash dilutes the liquid, but only small changes to the pH value are obtained (Figure 1).

Figure 1: Neutralizing effect of Plum pH Neutral and 0.9 % sodium chloride*



1: Prof. Dr. med. Schrage N.F., Facharzt für Augenheilkunde an der Augenklinik Merheim, Ostmerheimer Str. 200 in 51109 Köln, in Kooperation mit dem Aachener Centrum für Technologietransfer in der Ophthalmologie ACTO, Geschäftsstelle Karlsburgweg 9 in 52070 Aachen, Vortrag: Augenverätzungen. Was gibt es Neues? Regionalforum Arbeitsmedizin, Stuttgart 23.02.2006, Folien 44 ff.
2: Prof. Dr. med. Schrage N.F., Facharzt für Augenheilkunde an der Augenklinik Merheim, Ostmerheimer Str. 200 in 51109 Köln, in Kooperation mit dem Aachener Centrum für Technologietransfer in der Ophthalmologie ACTO, Geschäftsstelle Karlsburgweg 9 in 52070 Aachen, Vergleichsergebnisse von Plum pH Neutral an verätzten Augen, Spülversuche an mit NaOH verätzten Kaninchenaugen im Vergleich von Plum pH Neutral mit physiologischer Kochsalzlösung, 24.05.2005.



**200 ML
PLUM PH NEUTRAL
SINGLE CARTON**

Bottles with a 4.9% sterile phosphate buffer solution that quickly neutralizes acids and alkali. They can be used individually or as refills in eyewash stations and boxes.

Single carton with 200 ml Plum pH Neutral. The carton protects the bottle and is suitable for refilling stations and for mobile workplaces.

Content: 200 ml
Packaging: 10 x 200 ml
Article no: 4755

**200 ML
PLUM PH NEUTRAL**

The 200 ml bottle is easy to carry in the specially designed belt bag, in the first aid box, in the toolbox, etc. Delivered with dust cap.

Content: 200 ml
Packaging: 10 x 200 ml
Article no: 4752

**500 ML
PLUM PH NEUTRAL DUO**

The 500 ml bottle comes with the ergonomic DUO eyecup, which makes it possible to rinse both eyes at the same time.

Content: 500 ml
Packaging: 6 x 500 ml
Article no: 4801



**BELT BAG
FOR 200 ML BOTTLES**

Carry your eyewash in this practical bag for your belt, work clothes, etc.

The bag is easy to open and provides quick first aid in case of accidents. A safe and hygienic way to store and protect the bottle. Bottle is not included.

Content: -
Packaging: 10 units
Article no: 4692



**” ACCIDENTS INVOLVING
ACIDS AND ALKALI WILL
HAPPEN WHEN YOU
LEAST EXPECT THEM**

Plum pH Neutral DUO saved the eyesight of a SITA France employee. They handle plastic waste like plastic bottles with cleaning detergents. The liquid from a cleaning bottle went into both employee's eyes.

Thanks to Plum pH Neutral DUO, the harmful substance was neutralized, and the employee was already back at work the next day.



Plum Wound and Eyewash Sprays are portable, easy to use and saline-based sprays. They work at any angle and provide a balanced and soft flow rate. Ideal for the industrial sector, production industry, laboratories, transport of dangerous goods, first responders, law enforcement, firefighters, military and first rescue services.

Wound & Eyewash Spray

Quick and easy to use within seconds
- for skin and eyes!

Gentle wound and eyewash spray used for minor accidents like cuts and scrapes on the skin or rinsing dust and dirt from wounds and eyes. The wound and eyewash spray is available in a saline solution and a phosphate buffered saline solution. They offer safe rinsing with a quick and easy single hand operation. The dust cap protects the spray nozzle from contamination, the spray is suitable for multiple use and the can is fully recyclable.

A quick and easy solution that minimizes the risk of further infection or damage.



Cleanse the nozzle before use by briefly spraying in the air.



Place the nozzle minimum 10 cm from the wound or eye. Press the nozzle and spray. Avoid letting the nozzle touch the wound or the eye.





**50 ML / 250 ML
WOUND AND EYEWASH SPRAY**

A 0.9% sodium chloride solution that gives immediate care to minor accidents such as cuts and scrapes on the skin or dust and dirt in the eyes. Helps minimize the risk of further infection or damage.

Wound and Eyewash Spray

Content: 50 ml
Packaging: 24 x 50 ml
Article no: 45530
Languages: en, de, sv, da, fr, no
 fi, pt, it, nl, es, pl

Content: 50 ml
Packaging: 24 x 50 ml
Article no: 45531
Languages: en, tr, cs, sl, hu, sk,
 lt, et, lv, el, ru

Content: 250 ml
Packaging: 12 x 250 ml
Article no: 4554
Wall bracket: 4123



**200 ML
WOUND AND EYEWASH BUFFER SPRAY**

A 0.9% phosphate buffered sodium chloride solution for rinsing away dust and dirt from wounds and eyes. Wound healing is supported by this saline spray, which mechanically removes foreign bodies, reduces pathogens, as well as blood, clot and dead tissue and moistening the tissue in the process.

Wound and Eyewash Buffer Spray

Content: 200 ml
Packaging: 12 x 200 ml
Article no: 4556
Wall bracket: 4123

Plum Wound and Eyewash Spray is CE marked according to the Guidelines for Classification of Medical Device Directive (MDD 93/42/EEC)



**200 ML / 250 ML
WALL BRACKET**

Make your first aid visible by mounting the Wound and Eyewash Spray cans on the wall. This simple wall bracket will save you minutes if a minor eye injury or damages on the skin occur. Then you do not have to search for your Wound and Eyewash Spray but it will be close at hand.

Bracket for 200 ml and 250 ml spray

Content: 1
Packaging: 24 pcs
Article no: 4123

Product benefits

- Sterile until empty - allows repeated use
- Quick and easy single hand operation
- Safe rinsing
- Dust cap to protect the spray nozzle from contamination
- Easy to transport, classed as a non-flammable aerosol
- Up to 100% product emptying
- The can is fully recyclable
- Hypoallergenic
- Sterile solution without preservatives
- Shelf-life of 3 years
- Expiry date on the bottom of the can



Open eyewash stations

An open wall-mounted solution gives a good overview and quick and easy access to first aid

PROTECT YOUR COLLEAGUES, MAKE FIRST AID AVAILABLE

When an accident has happened, it is too late to look for help in a box or in a cupboard.

Systematize first aid with a wall-mounted station and make it easy to access when an accident happens. Choose between open wall stations and box solutions. The latter are particularly suitable for dusty work environments to protect the bottles and avoid adding more dust to an accident.

Make sure to have additional bottles at hand so that rinsing can continue on the way to further treatment.

EASY TO MOUNT

All wall stations are easy to mount and are delivered with pictograms that clearly shows how to use the products.



PLUM EYEWASH STATION

Wall station with 500 ml Plum Eyewash. Particularly suitable for small and mobile workplaces.



PLUM EYEWASH STATION DUO

Wall station with 1000 ml Plum Eyewash DUO. Particularly suitable for small and mobile workplaces where prolonged rinsing of both eyes may be required.

Packaging:	1
Article no:	4611
Dimensions:	
Eyewash plan	h: 29 x w: 22,8 cm
Mounting plate	h: 29 x w: 13,5 cm

Packaging:	1
Article no:	4802
Dimensions:	
Eyewash plan	h: 29 x w: 22,8 cm
Mounting plate	h: 29 x w: 13,5 cm



PLUM COMBI-STATION

Wall station with 200 ml Plum pH Neutral and 500 ml Plum Eyewash. Particularly suitable for places where there is a risk of getting both foreign bodies and acids and alkali in the eyes.

Packaging:	1
Article no:	4770
Dimensions:	
Eyewash plan	h: 29 x w: 22,8 cm
Mounting plate	h: 29 x w: 22,8 cm



PLUM COMBI-STATION DUO

Wall station with 500 ml Plum pH Neutral DUO and 1000 ml Plum Eyewash DUO. Particularly suitable for places where there is a risk of getting both foreign bodies and acids and alkali in the eyes and where there is a risk that both eyes will be affected.

Packaging:	1
Article no:	4803
Dimensions:	
Eyewash plan	h: 29 x w: 22,8 cm
Mounting plate	h: 29 x w: 22,8 cm



PLUM EYEWASH STATION

Wall station with 2 x 500 ml Plum Eyewash. Particularly suitable for environments and small workplaces that are not dirty, but where there is a risk of getting foreign bodies into the eyes.

Packaging:	1
Article no:	4694
Dimensions:	
Eyewash plan	h: 29 x w: 22,8 cm
Mounting plate	h: 29 x w: 22,8 cm



PLUM EYEWASH STATION MAXI

Wall station with 2 x 1000 ml Plum Eyewash. Particularly suitable for places where prolonged rinsing may be required.

Packaging:	1
Article no:	4708
Dimensions:	
Eyewash plan	h: 29 x w: 22,8 cm
Mounting plate	h: 29 x w: 22,8 cm

The green Plum Eyewash helps you remove dirt and foreign objects

The blue Plum pH Neutral provides immediate first aid in case of serious accidents involving acids and alkali

The soft bottles are easy to handle and rinse the eye with a soft and even flow

Clear pictograms and symbols show how to store and handle the bottle

The ergonomic soft plastic eyecup is shaped to fit the curve of the eye. Squeeze it together before putting it to your eye. The cup helps keep the eye open and provide the best possible flow of liquid and rinsing time



Quick and easy to open. A single turn breaks the seal and you are ready to rinse



Clear instructions for use

Economic and a long shelf-life of 3,5 years. The expiry date is clearly shown on the front of the label

With the ergonomic DUO eyecup, both eyes can be rinsed at the same time. The cup helps keep the eyes open and provide the best possible flow of liquid and rinsing time

Quickly accessible eyewash solutions

Quick access to eyewash is crucial for all workplaces with a risk of eye injuries, because seconds decide the impact of an accident

YOU ONLY HAVE TWO EYES – TAKE CARE OF THEM

Your eyesight is one of the most important senses. So, take care of your eyes. Accidents can happen no matter how careful you are. Stay prepared for accidents with the right first aid at the right place. An eyewash solution placed in high risk areas is the best preparation for eye injuries.

A WIDE RANGE OF UNIQUE RINSING SOLUTIONS

Our top priority is to help you get a safer place to work, and our concepts is developed in close collaboration with skilled eye specialists. This means that the products have been thought through down to the smallest detail. They are easy to use with clear pictograms that show you the correct use. The eyecup helps keep the eye open. The rinsing liquid hits the eye in the right place and with the right amount in a soft and even flow.



The products are CE marked according to the Guidelines for Classification of Medical Device Directive (MDD 93/42/EEC)

Plum iBox


– designed for your safety

Delivered with wall bracket and wall mounting set

Quick and efficient overview of eyewash bottles and expiry dates

Release from wall bracket and quickly bring to injured person

Quickly accessible eyewash solutions



Delivered with an eyewash plan

Strong and durable material (EPP)

Spacious design lets you grab bottles fast and easily

Easily opens - mounted on the wall or standing on the floor



Dust-proof and portable eyewash box



PLUM IBOX 1

Dust-proof and portable eyewash box that fits one eyewash bottle of your choice. Suitable for dusty environments and small workplaces, with a risk of eye injuries.

Delivered empty

Packaging:	1
Article no:	4900
Dimensions:	
Eyewash plan	h: 29 x w: 22,8 cm
iBox 1	h: 42,5 x w: 18,6 x d: 13,5 cm

Plum iBox balances key features in a modern design and is built in a strong and durable material. Designed to accommodate the need for urgent action in case of an emergency by combining the quick product overview from an open station with the dust-proof features from an eyewash box, Plum iBox suits a wide range of workplaces.

A clear see-through window lets you check the bottles' condition and keep track of expiry dates immediately. A spacious design lets you grab bottles fast and easily just as the box opens smoothly - mounted on the wall or standing on the floor. Plum iBox is delivered with wall bracket, eyewash plan and wall mounting set.

Install Plum iBox as close as possible to an exposed area or work station and according to the risk assesment. Release Plum iBox from the wall and quickly bring the eyewash to the injured person. Every second counts!



PLUM IBOX 1

Dust-proof and portable eyewash box with 1000 ml Plum Eyewash DUO. Suitable for dusty environments and small workplaces, with a risk of getting foreign bodies into the eyes.

Packaging:	1
Article no:	4901
Dimensions:	
Eyewash plan	h: 29 x w: 22,8 cm
iBox 1	h: 42,5 x w: 18,6 x d: 13,5 cm



PLUM IBOX 1

Dust-proof and portable eyewash box with 500 ml Plum Eyewash. Suitable for dusty environments and small workplaces, with a risk of getting foreign bodies into the eyes.

Packaging:	1
Article no:	4902
Dimensions:	
Eyewash plan	h: 29 x w: 22,8 cm
iBox 1	h: 42,5 x w: 18,6 x d: 13,5 cm



PLUM IBOX 2

Dust-proof and portable eyewash box that fits two eyewash bottles of your choice. Suitable for dusty environments and small workplaces, with a risk of eye injuries.

Delivered empty

Packaging:	1
Article no:	4920
Dimensions:	
Eyewash plan	h: 29 x w: 22,8 cm
iBox 2	h: 42,5 x w: 29 x d: 13,5 cm



PLUM IBOX 2

Dust-proof and portable eyewash box with 200 ml pH Neutral and 500 ml Plum Eyewash. Particularly suitable for places where there is a risk of getting both foreign bodies and acids and alkali in the eye.

Packaging: 1
Article no: 4921
Dimensions:
 Eyewash plan h: 29 x w: 22,8 cm
 iBox 2 h: 42,5 x w: 29 x d: 13,5 cm



PLUM IBOX 2

Dust-proof and portable eyewash box with 2 x 500 ml Plum Eyewash. Particularly suitable for dusty environments and small workplaces, with a risk of getting foreign bodies into the eyes.

Packaging: 1
Article no: 4922
Dimensions:
 Eyewash plan h: 29 x w: 22,8 cm
 iBox 2 h: 42,5 x w: 29 x d: 13,5 cm



PLUM IBOX 2

Dust-proof and portable eyewash box with 500 ml Plum pH Neutral DUO and 1000 ml Plum Eyewash DUO. Particularly suitable for places with a risk of getting both foreign bodies and acids and alkali in the eyes and where there is a risk that both eyes will be affected.

Packaging: 1
Article no: 4923
Dimensions:
 Eyewash plan h: 29 x w: 22,8 cm
 iBox 2 h: 42,5 x w: 29 x d: 13,5 cm



PLUM IBOX 2

Dust-proof and portable eyewash box with 2 x 1000 ml Plum Eyewash DUO. Particularly suitable for places where prolonged rinsing may be required.

Packaging: 1
Article no: 4924
Dimensions:
 Eyewash plan h: 29 x w: 22,8 cm
 iBox 2 h: 42,5 x w: 29 x d: 13,5 cm



PLUM IBOX 3

Dust-proof and portable eyewash box that fits three eyewash bottles and one spray. Suitable for dusty environments and small workplaces, with a risk of eye injuries.

Delivered empty

Packaging: 1
Article no: 4930
Dimensions:
 Eyewash plan h: 29 x w: 22,8 cm
 iBox 3 h: 47,5 x w: 37,6 x d: 13,5 cm



PLUM IBOX 3

Dust-proof and portable eyewash box with 200 ml pH Neutral, 500 ml Plum Eyewash, 1000 ml Plum Eyewash DUO and 250 ml Wound and Eyewash Spray. Particularly suitable for places where there is a risk of getting both foreign bodies and acids and alkali in the eye.

Packaging: 1
Article no: 4931
Dimensions:
 Eyewash plan h: 29 x w: 22,8 cm
 iBox 3 h: 47,5 x w: 37,6 x d: 13,5 cm





Dust-proof eyewash box

A closed station protects your first aid equipment from dirt and dust and keeps the products ready for use

PLUM COMBI-BOX

Dust-proof eyewash box with 200 ml pH Neutral and 500 ml Plum Eyewash. Suitable for places where there is a risk of getting both foreign bodies and acids and alkali in the eye.

Packaging: 1
Article no: 4789
Dimensions:
 Eyewash plan h: 29 x w: 22,8 cm
 Eyewash box h: 28 x w: 23 x d: 11 cm

REDUCE IMPACT WITH COMBINED RISING

Combined rinsing is recommended in case of accidents involving acids and alkali – first with Plum pH Neutral, then with Plum Eyewash. A combi-box holds both products and gives quick and easy access. Start by neutralizing the acid or alkali with Plum pH Neutral, continue rinsing with Plum Eyewash until you reach a doctor or the emergency room.

If you work mostly with liquid chemicals, we recommend Plum pH Neutral DUO. That way you can rinse both eyes at the same time in case of an emergency.

We recommend the open stations for environments without dust and dirt, such as laboratories. The closed stations protect the products in dusty environments. All stations are easy to mount and come with a mirror and pictograms clearly indicating how to use the eyewash bottles.



PLUM EYEWASH BOX

Dust-proof eyewash box with 2 x 500 ml Plum Eyewash. Particularly suitable for dusty environments and small workplaces, with a risk of getting foreign bodies into the eyes.

Packaging: 1
Article no: 4650
Dimensions:
 Eyewash plan h: 29 x w: 22,8 cm
 Eyewash box h: 28 x w: 23 x d: 11 cm



PLUM COMBI-BOX DUO

Dust-proof eyewash box with 500 ml Plum pH Neutral DUO and 500 ml Plum Eyewash DUO. Particularly suitable for places with a risk of getting both foreign bodies and acids and alkali in the eyes and where there is a risk that both eyes will be affected.

Packaging: 1
Article no: 4862
Dimensions:
 Eyewash plan h: 29 x w: 22,8 cm
 Eyewash box h: 28 x w: 23 x d: 11 cm



PLUM COMBI-BOX DUO MAXI

Dust-proof eyewash box with 500 ml Plum pH Neutral DUO and 1000 ml Plum Eyewash DUO. Particularly suitable for places with a risk of getting both foreign bodies and acids and alkali in the eyes and where there is a risk that both eyes will be affected.

Packaging: 1
Article no: 4810
Dimensions:
 Eyewash plan h: 29 x w: 22,8 cm
 Eyewash box h: 35,7 x w: 25,2 x d: 11,6 cm

The soft close system makes it easy to open and close the eyewash box and protects it from dust and dirt from the surroundings

Easy to mount close to the working station

The moulded holders fit all Plum Eyewash and Plum pH Neutral bottles - both with regular and DUO eyecups - and keep them securely stored

Made of polystyrene - a light, isolating and environmentally friendly material. It is produced with less plastic compared to a 100% plastic construction.

When replaced the eyewash box is suitable for burning

Appointed "Worldstar Winner" for being an innovative product because the packaging is an integrated and permanent part of the product



ENSURE THE RIGHT TEMPERATURE
 Whatever the reason, most often quick help is needed when you get something in the eye. Therefore, preparedness must be in place and be easily accessible at the workplace – also under cold weather conditions.

PLUM EYEWASH BOX DUO MAXI

Dust-proof eyewash box with 2 x 1000 ml Plum Eyewash DUO. Particularly suitable for places where prolonged rinsing may be required.

PLUM EYEWASH PLUG & HEAT BOX, SMALL

Dust-proof eyewash box with a heating element. Ensures the eyewash liquid has a comfortable shower temperature. Suitable for 2 bottles of 200 ml and/or 500 ml Plum Eyewash or pH Neutral. No content.

Packaging:	1
Article no:	4816
Dimensions:	
Eyewash plan	h: 29 x w: 22,8 cm
Eyewash box	h: 35,7 x w: 25,2 x d: 11,6 cm

Packaging:	1
Article no:	4672
Dimensions:	
Eyewash plan	h: 29 x w: 22,8 cm
Eyewash box	h: 28 x w: 23 x d: 11 cm



For fast and effective decontamination of bodies and surfaces in contact with corrosive and toxic hazardous substances such as strong acids and bases, oxidants, formaldehyde, phenols, organic solvents, hydrofluoric acid, and fluorides.

Plum BodyNeutrAll

Mobile body shower with a highly effective decontamination solution for skin, surfaces and clothing.

Plum BodyNeutrAll contains ascorbic acid, an amino acid buffer, and calcium ions as active substances. Glycerin, polyethylene and glycol are solubilizers for organic chemicals.

PLUM BODYNEUTRALL:

- ✓ reduces strong oxidants (such as chlorine, hypochlorite, chlorine dioxide, chloramine, bromine, iodine, fluorine, hydrogen peroxide, peracetic acid) to compounds that are no longer active.
- ✓ binds aldehydes into imines, phenols into amines, and heavy metals into ammine complexes.
- ✓ binds hydrofluoric acid and fluorides into the sparingly soluble calcium fluoride.
- ✓ dissolves organic solvents and chemicals.

PERMITS / APPROVALS OF THE DECONTAMINATION SOLUTION:

- Certification as a class I medical device. according to the requirements of the EU Regulation EU 2017/745.
- Officially registered by the Federal Office for Drugs and Medical Devices under the number DE/CA33/ILM/2020/1. The registration is valid for the entire EU and may be marked with the CE mark accordingly.
- Conforms to DIN 14555-12 for hazardous goods equipment trolley.



PLUM BODYNEUTRALL 5 L

5 L empty sprayer which comes with a canister of mixed liquid to fill in the sprayer. This is suitable for larger areas up to 50 m² and covers several people.

Packaging: 1
Article no: 4735
Dimensions:
 Sprayer h: 57 (with lance: 80) x w: 22 cm
 Ø18 cm



PLUM BODYNEUTRALL 600 ML

2 bottles of 300 ml to be mixed and used with a sprayhead. Easy to use on smaller areas up to 6 m² or up to 3 persons.

Packaging: 4
Article no: 4736
Dimensions:
 2 bottles with sprayhead h: 29 x w: 11 x d: 6 cm



PLUM BODYNEUTRALL 100 ML

Mixed, pocket-sized bottle to bring along and it covers approx. 1 m² or 1 arm's length.

Packaging: 6
Article no: 4737
Dimensions:
 h: 13 x w: 6.5 x d: 3.1 cm



C-GEL 30 G

Prevents the effect of toxic fluorides on the skin in case of contamination with hydrofluoric acid. Simplified application in an emergency - can also be used with gloves - due to the pressure dispenser.

Packaging: 8
Article no: 4738



For further information, please scan the QR code

Quick

FIX

QUICKFIX IS A UNIQUE DISPENSER:

- Immediate and visible access to plasters and eyewash ampoules
- Each plaster is ready for use with a pull from the dispenser
- Each plaster is packed separately
- Only one hand is needed to apply
- Easily refilled as required
- Easy to see when it is empty
- Quickly stops minor bleeding without staining the other plasters
- Small dimensions - easy to find room for the QuickFix dispenser
- Includes a shelf for storing small, relevant objects
- Can be locked so that it can only be opened with a key



Your quick help for minor everyday accidents

Mount the plaster dispenser on the wall and save time.
Put on the plaster - and back to work!

QUICK AND HYGIENIC HELP WITH MINOR INJURIES

Stop wasting time on minor injuries. Take a plaster from the dispenser, put it on, and go back to work! QuickFix solves the problem easily, hygienically and effectively.

Lots of minor accidents happen every day at work where a plaster is all that is needed. Choose the right plaster for your work situation. A suitable solution is available whether you need plasters for common work situations, for workplaces with a lot of humidity or for work within the food sector.



Take the plaster from the dispenser.
It is ready for use...



... and can be put on
using one hand



How to use QuickFix
(animation video)



Flexible system for all types of work

QuickFix provides a structured approach
to first aid of minor injuries

CHOOSE THE DISPENSER THAT BEST SUITS YOUR NEEDS

Our range of QuickFix plaster dispensers contain two types of dispensers in different designs. They are developed for different types of workplaces and they can be filled with plasters and ampoules of eyewash in a variety of mixes to meet the needs of every work situation.




QUICKFIX DISPENSER

The classic green QuickFix dispenser contains two refills of QuickFix plasters and/or QuickRinse eyewash ampoules. The green cap protects the content of the dispenser from dust and dirt from the workplace. It is designed with an inside shelf for smaller objects and it can be secured with lock and key.

The dispenser can be used in any location with risk of injuries and scratches. The dispenser is designed for mounting on any stationary workplace.

Dimensions: h: 135 x w: 230 x d: 32 mm



QUICKFIX UNO DISPENSER

For workplaces that are not challenged with dust and dirt, we have developed the QuickFix UNO plaster dispenser. It is a small and modern dispenser in a unique design that combines workplace safety with modern style.

The dispenser contains 1 refill of QuickFix plasters or QuickRinse eyewash ampoules.

QuickFix UNO is available with a glossy white or a transparent cover.

Dimensions: h: 85 x w: 130 x d: 35 mm

Choose your QuickFix plaster dispenser

- to meet your specific needs



QUICKFIX WATER RESISTANT
 Dispenser with 90 water resistant plasters suitable for environments with a high level of humidity.

Packaging: 10
Article no: 5501
Refill article no: 5511



QUICKFIX ELASTIC & WATER RESISTANT
 Dispenser with 45 elastic plasters and 45 water resistant plasters.

Packaging: 10
Article no: 5507
Refill article no: 5512 & 5511



QUICKFIX ELASTIC
 Dispenser with 90 elastic textile plasters that adjust to the movement of the skin.

Packaging: 10
Article no: 5502
Refill article no: 5512



QUICKFIX ELASTIC LONG
 Dispenser with 60 long elastic plasters that can be wrapped around the finger several times.

Packaging: 10
Article no: 5528
Refill article no: 5508



QUICKFIX DETECTABLE
 Dispenser with 90 blue plasters with built-in metal surface which can be traced by a metal detector.

Packaging: 10
Article no: 5503
Refill article no: 5513



QUICKFIX DETECTABLE LONG
 Dispenser with 60 long blue plasters with built-in metal surface which can be traced by a metal detector.

Packaging: 10
Article no: 5529
Refill article no: 5509



QUICKFIX ALOE VERA
 Dispenser with 90 aloe vera plasters, which soothes and moisturizes.

Packaging: 10
Article no: 5506
Refill article no: 5514



QUICKFIX&RINSE
 Dispenser with 45 elastic plasters and 5 ampoules with 20 ml 0.9% sterile sodium chloride.

Packaging: 10
Article no: 4631
Refill article no: 5512 & 5160



QUICKFIX ALU
 Dispenser with 90 alu plasters with anti-adhering effect due to the aluminum layer.

Packaging: 10
Article no: 5505
Refill article no: 5515



QUICKRINSE
 Dispenser with 10 eyewash ampoules with 20 ml 0.9% sterile sodium chloride.

Packaging: 10
Article no: 4632
Refill article no: 5160



QUICKFIX BLOOD STOPPER
 Dispenser with 90 blood stopper plasters with hemostatic effect.

Packaging: 10
Article no: 5510
Refill article no: 5516



QUICKFIX EMPTY DISPENSER
 Dispenser without refills. Fill with any QuickFix or QuickRinse refills of your choice.

Packaging: 10
Article no: 5500



**QUICKFIX UNO ELASTIC
DISPENSER, TRANSPARENT**

Transparent dispenser with 45 elastic and comfortable textile plasters that adjust to the movement of the skin. The plasters are skin-friendly and allow the skin to breathe.

Packaging: 10
Article no: 5531

**QUICKFIX UNO ELASTIC
DISPENSER, WHITE**

Glossy white dispenser with 45 elastic and comfortable textile plasters that adjust to the movement of the skin. The plasters are skin-friendly and allow the skin to breathe.

Packaging: 10
Article no: 5532



**QUICKFIX UNO DISPENSER
WITHOUT REFILLS, WHITE**

Fill the dispenser with the plaster of your own choosing. All the QuickFix plaster refills can be used according to your specific needs. The dispenser can also be used with a QuickRinse refill.

Packaging: 10
Article no: 5534



**QUICKFIX UNO DISPENSER
WITHOUT REFILLS, TRANSPARENT**

Packaging: 10
Article no: 5533



QUICKFIX MINI

Small and practical pack with 30 elastic plasters. Easy to carry along anywhere you may need a plaster. Each plaster is packed separately and is ready for use just by pulling it from the pack.

Packaging: 30 x 30
Article no: 5504
Plaster dimensions: l: 72 x w: 25 mm



QuickFix plaster refills

Full overview of refill solutions



QUICKFIX ELASTIC REFILL

45 elastic and comfortable textile plasters that adjust to the movements of the skin. The plasters are skin-friendly and allow the skin to breathe.



QUICKFIX DETECTABLE REFILL

45 blue conspicuous plasters with a built-in metal foil that can be traced by a metal detector. Particularly suitable for the food industry. The plasters are made from a breathable material.



QUICKFIX WATER RESISTANT REFILL

45 water resistant plasters that are particularly suitable for working environments with a high level of humidity. The plasters are made from a protective PE material.

REFILL

Display carton: 6 x 45
Packaging: 8 x 6 x 45
Article no: 5512
Dimensions: l: 72 x w: 25 mm

DISPENSER

Article no: 5502

REFILL

Display carton: 6 x 45
Packaging: 8 x 6 x 45
Article no: 5513
Dimensions: l: 72 x w: 25 mm

DISPENSER

Article no: 5503

REFILL

Display carton: 6 x 45
Packaging: 8 x 6 x 45
Article no: 5511
Dimensions: l: 72 x w: 25 mm

DISPENSER

Article no: 5501



QUICKFIX LONG ELASTIC AND DETECTABLE

30 plasters with extra length that can be wrapped around the finger several times. Choose between elastic textile plasters or blue detectable plasters.



ELASTIC LONG REFILL

Display carton: 6 x 30
Packaging: 8 x 6 x 30
Article no: 5508
Dimensions: l: 120 x w: 20 mm

DISPENSER

Article no: 5528

DETECTABLE LONG REFILL

Display carton: 6 x 30
Packaging: 8 x 6 x 30
Article no: 5509
Dimensions: l: 120 x w: 20 mm

DISPENSER

Article no: 5529



QUICKFIX ALOE VERA REFILL

The aloe vera contained in the pad has well known soothing, moisturizing and calming properties. It creates the best conditions for the wound environment helping the natural healing process. Prevents the formation of scabs.

REFILL

Display carton: 6 x 45
Packaging: 8 x 6 x 45
Article no: 5514
Dimensions: l: 72 x w: 25 mm

DISPENSER

Article no: 5506



QUICKFIX ALU REFILL

Plaster with micronized aluminum pad. The anti-adhering layer in aluminum and PE has a higher anti-adhering performance compared to a standard PE layer. The plaster is air-permeable and comfortable.

REFILL

Display carton: 6 x 45
Packaging: 8 x 6 x 45
Article no: 5515
Dimensions: l: 72 x w: 25 mm

DISPENSER

Article no: 5505



QUICKFIX BLOOD STOPPER REFILL

Plasters with hemostatic effect. The presence of alginate fibres between the nonwoven fibres of the wound pad creates an artificial clot (gel) inside the pad, which blocks the bleeding in 1-2 minutes and maintains optimal moist conditions.

REFILL

Display carton: 6 x 45
Packaging: 8 x 6 x 45
Article no: 5516
Dimensions: l: 72 x w: 25 mm

DISPENSER

Article no: 5510



5 X 20 ML QUICKRINSE REFILL

5 eyewash ampoules with 20 ml 0.9% sterile sodium chloride. For use individually or for refilling QuickSafe®, QuickFix&Rinse and QuickRinse dispensers.

Content: 5 x 20 ml
Packaging: 4 x 5 x 20 ml
Article no: 5160



QuickSafe® gives you quick, easy and safe first aid

A flexible solution that gives a quick overview of the available first aid - and makes it easy to always have the right products in the right place

CUSTOMIZED FIRST AID STATIONS

Fill your first aid station according to your needs. Eyewash for rinsing foreign bodies from the eye. Plum pH Neutral for accidents with acids or alkali. Plasters that are suitable for your work situation. Burn gel that soothes minor burns. Eyewash ampoules for minor accidents and irritated eyes. Wipes for cleansing wounds. Compress dressings for more serious accidents.

All of it behind a transparent cover that will give you a perfect overview of what you need.

” GET THE MOST IMPORTANT ITEMS IN THE RIGHT PLACE

Spæncom A/S is a leading supplier of concrete elements. The company has chosen Plum Quicksafe® in order to increase safety for their more than 500 employees. “We have chosen Plum Quicksafe®, because the station can be organized to suit the individual work area. In the areas where concrete elements are acid-washed, we can, for example, place Plum Eyewash and Plum pH Neutral DUO close to the work situation. QuickSafe® contains the most necessary items, it can be placed so that first aid is close at hand, and it is very practical with the transparent cover.”

Jany Johannsen,
Environment and Safety Coordinator





QuickSafe® EMPTY

YOU MAKE THE DECISIONS

If you want full freedom to organize your first aid according to your requirements, you can start with an empty station and fill it to suit your specific situation.



BUY AN EMPTY STATION AND FILL IT TO MEET YOUR NEEDS

If you have requirements that we do not meet with the 4 different QuickSafe® solutions presented on the following page, you can get a 100% customized solution. Start with a wall station and fill it with eye-wash bottles and other first aid products.


Packaging: 1
Article no: 5173
Dimensions: h: 430 x w: 253 x d: 92 mm

Help is at hand

Easy access to first aid is the key to increased safety

QuickSafe® collects the most important first aid products in one place. Position the station close to where the products are most likely to be used. The station does not take up much more room than an A3 sheet of paper, which means you can place your QuickSafe® exactly where it is needed.

The transparent cover gives you an overview and provides quick and easy access to first aid. You can easily see when the station needs a refill, and the dust-proof cover together with the separately packed products make QuickSafe® both a hygienic and an elegant solution.



The eyewash plan is included in all QuickSafe® stations



QuickSafe® BASIC

A GOOD STARTING POINT

With a QuickSafe® Basic you get the most essential products and can choose the contents for the top shelves. It gives you full flexibility.



DO YOU HAVE SPECIAL REQUIREMENTS REGARDING YOUR FIRST AID STATION?

The backbone of the entire QuickSafe® series is Plum's acclaimed eyewash. With QuickSafe® Basic you will get it in a combination with Plum pH Neutral, eyewash ampoules and plasters. The rest is up to you. Organize the shelf in the station to suit your company's requirements.

CONTENTS:	Article no:
2 x 500 ml Plum Eyewash	4604
1 x 200 ml Plum pH Neutral	4752
5 x 20 ml QuickRinse eyewash ampoules	5160
1 x QuickFix Elastic refill with 45 plasters	5512

Packaging:	1
Article no:	5170
Dimensions:	h: 430 x w: 253 x d: 92 mm



QuickSafe® COMPLETE

THE COMPLETE SOLUTION

Eyewash including Plum pH Neutral and eyewash ampoules, pressure dressing, wound cleansing wipes, plasters and burn gel to soothe minor burns.



WIDE FIRST AID COVER

With QuickSafe® Complete you will get a wide range of first aid products. There are products for most critical situations and for many of the everyday minor accidents.

CONTENTS:	Article no:
18 x QuickCool burn gel	5150
1 x QuickStop wound compress dressing 1 x 17 x 17 cm and 2 x 12 x 7.5 cm	5152
20 x QuickClean wound cleanser	5151
2 x 500 ml Plum Eyewash	4604
1 x 200 ml Plum pH Neutral	4752
5 x 20 ml QuickRinse eyewash ampoules	5160
1 x QuickFix Elastic refill with 45 plasters	5512

Packaging:	1
Article no:	5174
Dimensions:	h: 430 x w: 253 x d: 92 mm



QuickSafe® FOOD INDUSTRY

FOR THE SAFETY OF YOU AND YOUR CUSTOMERS

We have filled a station with the relevant content to suit the requirements of the food industry.



QUICK AND SAFE HELP WITH BUILT-IN TRACEABILITY

A first aid station targeting the food industry, kitchens, canteens and restaurants with, among other things, traceable plasters, wound cleansing, burn gel and wound compress dressings.

CONTENTS:	Article no:
18 x QuickCool burn gel	5150
1 x QuickStop wound pressure dressing	5152
1 x 17 x 17 cm and 2 x 12 x 7.5 cm	
20 x QuickClean wound cleanser	5151
2 x 500 ml Plum Eyewash	4604
1 x 200 ml Plum pH Neutral	4752
1 x QuickFix Detectable refill with 45 plasters	5513
1 x QuickFix Detectable Long refill with 30 plasters	5509

Packaging:	1
Article no:	5175
Dimensions:	h: 430 x w: 253 x d: 92 mm



QuickSafe® CHEMICAL INDUSTRY

STRONG SOLUTION FOR TOUGH ENVIRONMENTS

Accidents with acids/alkali can have serious and life-long consequences for the individual if you compromise on accessibility.



WHEN ONLY SECONDS STAND BETWEEN A MINOR AND A MAJOR ACCIDENT

If there is a risk of accidents involving acids or alkali, safety in the form of effective and mobile eyewash is crucial. It contains, among other things, Plum pH Neutral DUO that will rinse both eyes at the same time.

CONTENTS:	Article no:
1 x 500 ml Plum Eyewash DUO	4861
1 x 500 ml Plum pH Neutral DUO	4801
1 x 200 ml Plum pH Neutral	4752
5 x 20 ml QuickRinse eyewash ampoules	5160
1 x QuickFix Elastic refill with 45 plasters	5512

Packaging:	1
Article no:	5171
Dimensions:	h: 430 x w: 253 x d: 92 mm

Plum Pull1Aid®

Quick, easy and safe first aid

Plum Pull1Aid® is a wall mounted first aid station that fits most workplaces. A flexible solution that gives a quick overview of the available first aid - and makes it easy to always have the right products in the right place.

Pull1Aid® allows you to combine first aid and wound care products that suit the needs at your workplace. The flexible combinations provide a tailored solution and avoid excess products.



Plum Pull1Aid® Stations



CONTENTS:

1 x QuickCool burn gel	ART. NO.:	5150
1 x QuickClean wound cleanser		5151
1 x QuickStop wound pressure dressing		5152
2 x Blood Stopper wound pressure dressing		5154
1 x Blood Stopper Mini wound pressure dressing		5153
1 x Burn Gel dressing for burns		5156
1 x Protect Kit respiratory aid		5155
1 x 200 ml Wound- and Eyewash Spray		4556
1 x QuickFix Elastic		5512
1 x QuickRinse		5160

Packaging:	1
Article no:	5181
Dimensions:	h: 430 x w: 253 x d: 92 mm

CONTENTS:

1 x QuickCool burn gel	ART. NO.:	5150
1 x QuickClean wound cleanser		5151
2 x Blood Stopper wound pressure dressing		5154
1 x Blood Stopper Mini wound pressure dressing		5153
1 x Burn Gel dressing for burns		5156
1 x Protect Kit respiratory aid		5155
1 x 200 ml Wound- and Eyewash Spray		4556
1 x QuickFix Elastic		5512
1 x QuickFix Water Resistant		5511

Packaging:	1
Article no:	5180
Dimensions:	h: 430 x w: 253 x d: 92 mm

Plum Pull1Aid® Refills

BLOOD STOOPER MINI

A smaller, sterile compress dressing with four functions. Use as pressure pad for heavy bleeding, protection for superficial wounds, burn dressing with wound-protecting surface or as temporary support.



Art. no. 5153

BLOOD STOPPER

A sterile compress dressing with four functions. Use as pressure pad for heavy bleeding, protection for superficial wounds, burn dressing with wound-protecting surface or as temporary support.



Art. no. 5154

PROTECT KIT

Protecting aid for mouth-to-mouth resuscitation situations. Limits contact between the first aid provider and the person on harm by a breathing mask is with a one-way-valve.



Art. no. 5155

BURN GEL

Gel dressing that rapidly cools and effectively relief pain for burns. The cooling effect reduces pain and swelling and prevents the burn from going deeper into the skin.



Art. no. 5156

QUICKCOOL

Burn gel for first aid in case of minor burns and scalding. Cooling and pain-relieving effect. Contains tea tree oil. Can be used individually and for refilling the first aid stations.



Art. no. 5150

QUICKCLEAN

Wound cleansing wipes for cleaning minor wounds and cuts. Each wipe is packed separately and contains sterile water. Can be used individually and for refilling first aid stations.



Art. no. 5151

QUICKSTOP

Wound compress dressing set for minor as well as major bleeding. Sterile dressing with integrated compression. Elastic, non-adhesive and skin-friendly. For general use as pressure, support and protection dressing.



Art. no. 5152

WOUND AND EYEWASH SPRAY

A 0.9% phosphate buffered sodium chloride solution, that removes foreign bodies, blood, clot and dead tissue and gives immediate care to minor eye accidents or cuts and scrapes on the skin. 200 ml.



Art. no. 4556

QUICKFIX ELASTIC

Elastic and comfortable textile plasters that adjust to the movements of the skin. The plasters are skin-friendly and allow the skin to breathe.



Art. no. 5512
Art. no. long 5508

QUICKFIX DETECTABLE

Blue conspicuous plasters with a built-in metal foil that can be traced by a metal detector. Made from a breathable material.



Art. no. 5513
Art. no. long 5509

QUICKFIX WATER RESISTANT

Water resistant plasters that are suitable for working environments with a high level of humidity. The plasters are made from a protective PE material.



Art. no. 5511

QUICKFIX ALOE VERA

Aloe vera contained plasters with soothing, moisturizing and calming effect. Creates the best conditions for natural wound healing and prevents the formation of scabs.



Art. no. 5514

QUICKFIX ALU

Air-permeable and comfortable plaster with micronized aluminum pad. Anti-adhering layer in aluminum and PE with high anti-adhering performance.



Art. no. 5515

QUICKFIX BLOOD STOPPER

Hemostatic effect with alginate fibres between the nonwoven fibres of the wound pad. Creates artificial clot inside the pad, which blocks bleeding in 1-2 minutes and maintains optimal moist conditions.



Art. no. 5516

QUICKRINSE

5 eyewash ampoules with 20 ml 0.9% sterile sodium chloride. For use individually or for refilling Pull1Aid®, QuickSafe®, QuickFix&Rinse and QuickRinse dispensers.



Art. no. 5160

Wound care

Practical solutions
– to stop small accidents
before they stop you

Even small cuts and
wounds can get serious
if not tended to in time.
With the right wound
care at your disposal,
you can easily cleanse
and dress a small wound
- and quickly get back to
work.



FIRST AID BAG - BASIC

Practical first aid bag for your belt and workwear or any workplace requiring a portable solution for first aid. Easy to open and provides quick and hygienic first aid on the go.

BAG INCLUDES:

- 1 x QuickFix Water Resistant
- 1 x QuickFix Elastic
- 1 x Plum Eyewash 200 ml
- 1 x QuickCool
- 1 x QuickClean
- 1 x QuickStop

Packaging: 4
Article no: 4960



FIRST AID BAG - INDUSTRIAL

Practical first aid bag for your belt and workwear or any workplace requiring a portable solution for first aid. Easy to open and provides quick and hygienic first aid on the go.

BAG INCLUDES:

- 2 x QuickFix Elastic
- 1 x pH Neutral 200 ml
- 1 x Wound and eyewash spray 50 ml
- 1 x QuickCool
- 1 x QuickClean
- 1 x QuickStop

Packaging: 4
Article no: 4961



QUICKSTOP

Wound compress dressing set for effective first aid in the event of minor as well as major bleeding. Sterile dressing with integrated compression. Elastic, non-adhesive and skin-friendly. QuickStop for general use as pressure, support and protection dressing.

SET INCLUDES:

- 1 x wound compress dressing, 17 x 17 cm
- 2 x wound compress dressings, 12 x 7.5 cm

Packaging: 6 x 3
Article no: 5152



**50 ML / 250 ML
WOUND AND EYEWASH SPRAY**

A 0.9% sodium chloride solution that gives immediate care to minor eye accidents or cuts and scrapes on the skin.

Content: 50 ml
Packaging: 24 x 50 ml
Article no: 45530
Languages: en, de, se, da, fr, no, fi, pt, it, nl, es, pl

Content: 50 ml
Packaging: 24 x 50 ml
Article no: 45531
Languages: en, tr, cs, si, hu, sk, lt, et, lv, el, ru

Content: 250 ml
Packaging: 12 x 250 ml
Article no: 4554



**200 ML
WOUND AND EYEWASH BUFFER SPRAY**

A 0.9% phosphate buffered sodium chloride solution, that removes foreign bodies, blood, clot and dead tissue.

Content: 200 ml
Packaging: 12 x 200 ml
Article no: 4556

Plum Wound and Eyewash Spray is CE marked according to the Guidelines for Classification of Medical Device Directive (MDD 93/42/EEC)



ALOE VERA COOL GEL

Highly concentrated organic Aloe Vera gel that cools, hydrates, and softens the skin.

Aloe Vera is a purely natural product, has a regenerating effect on the skin and is suitable for all skin types.

Apply the soothing gel abundantly on insect bites or exposed areas such as sunburned skin.

Content: 75 ml
Packaging: 8 x 75 ml
Article no: 5570



PLUM NANO SAFE

For effective removal of nanoscale particles from the skin.

Content: 50 ml
Article no: 5576



QUICKCOOL

Burn gel for first aid in case of minor burns and scalding. Cooling and pain-relieving effect. Contains tea tree oil. Can be used individually and for refilling the QuickSafe® stations. 18 sachets per box.

Packaging: 6 x 18
Article no: 5150



QUICKCLEAN

Wound cleansing wipes for cleaning minor wounds and cuts. Each wipe is packed separately and contains sterile water. Can be used individually and for refilling the QuickSafe® stations. 20 sachets per box.

Packaging: 6 x 20
Article no: 5151



QUICKCLEAN WOUND CLEANSING

Box with 40 wound cleansing wipes for cleaning minor wounds and cuts. Each wipe is packed separately and contains sterile water. Each wipe measures 130 x 200 mm.

Packaging: 11 x 40
Article no: 5551
Wall bracket: 955904





Plum Burn Out

Fire Suppression Spray

Plum Burn Out is a fire suppression spray and a practical supplement to a regular fire extinguisher. The spray bottle is easy to bring along and contains an environmentally neutral liquid. It is highly effective, has fantastic reignition protection and it is safe to use.

- ✓ Environmentally neutral liquid
- ✓ Safe to use
- ✓ Mount on the wall in wall bracket
- ✓ Suitable for wood, paper and fabrics
- ✓ Empties in 20 seconds
- ✓ Scope 2 meters



PLUM BURN OUT


Content	250 ml
Packaging	12
Art. No.:	5590
Wall bracket:	4123

From hand care to eyewash

Plum Safety ApS departs from the Danish company Plum A/S – dating back to 1860

Plum A/S has been providing efficient and user-friendly hand hygiene solutions to corporate clients for more than 50 years – and for more than 20 years we have served the market with first-aid solutions that make your workplace a safer place to be.

Today we serve a wide range of users in the industry, handicraft and institutional areas. We offer first aid products such as eyewash, first aid stations, plasters, and burn gel.



Plum Safety ApS

Mandelalleen 1
DK-5610 Assens

T: +45 6916 9600
info@plum.eu

www.plum.eu

