



SAFETY AT WORK

# Plum Pull1Aid®



info@plum.eu  
+45 6916 9600

plum  
Safety



# Plum Pull1Aid®

Quick, easy and safe first aid

---

Easy access to first aid is the key to increase safety at work

Plum Pull1Aid® is a wall mounted first aid station that fits most workplaces. A flexible solution that gives a quick overview of the available first aid - and makes it easy to always have the right products in the right place.

Plum Pull1Aid® allows you to combine first aid and wound care products that suit the needs at your workplace. The flexible combinations provide a tailored solution and avoid excess products.

## PRODUCT FACTS & BENEFITS

- ✓ Clear cover - hygienic and protects the products from dust
- ✓ Practical and functional station
- ✓ High quality products
- ✓ Simple to use with clear instructions
- ✓ Easy and quick to refill
- ✓ Immediate and visible access to plasters, ready to use with a pull
- ✓ Wound care products are quick to pull out of the stations

Plum Pull1Aid® is clear, fast and simple first aid so you can act immediately if an accident happens. Expand your Plum Pull1Aid® station for increased safety for wounds and eyes by mounting a QuickSafe® station next to it.

# Plum Pull1Aid® Stations



## CONTENTS:

1 x Plum QuickCool burn gel	ART. NO.:	5150
1 x Plum QuickClean wound cleanser		5151
1 x Plum QuickStop wound dressing		5152
2 x Plum Blood Stopper wound dressing		5154
1 x Plum Blood Stopper Mini wound dressing		5153
1 x Plum BurnGel dressing for burns		5156
1 x Plum Protect Kit respiratory aid		5155
1 x 200 ml Plum Wound and Eyewash Spray		4556
1 x Plum QuickFix Elastic refill		5512
1 x Plum Eyewash 0.9% Sodium Chloride QuickRinse		5160

**Packaging:** 1  
**Article no:** 5181  
**Dimensions:** h: 430 x w: 253 x d: 92 mm

## CONTENTS:

1 x Plum QuickCool burn gel	ART. NO.:	5150
1 x Plum QuickClean wound cleanser		5151
2 x Plum Blood Stopper wound dressing		5154
1 x Plum Blood Stopper Mini wound dressing		5153
1 x Plum BurnGel dressing for burns		5156
1 x Plum Protect Kit respiratory aid		5155
1 x 200 ml Plum Wound and Eyewash Spray		4556
1 x Plum QuickFix Elastic refill		5512
1 x Plum QuickFix Water Resistant refill		5511

**Packaging:** 1  
**Article no:** 5180  
**Dimensions:** h: 430 x w: 253 x d: 92 mm

# Plum Pull1Aid® Refills

## PLUM BLOOD STOOPER MINI

A smaller, sterile compress dressing with four functions. Use as pressure pad for heavy bleeding, protection for superficial wounds, burn dressing with wound-protecting surface or as temporary support.



Art. no. 5153

## PLUM BLOOD STOPPER

A sterile compress dressing with four functions. Use as pressure pad for heavy bleeding, protection for superficial wounds, burn dressing with wound-protecting surface or as temporary support.



Art. no. 5154

## PLUM PROTECT KIT

Protecting aid for mouth-to-mouth resuscitation situations. Limits contact between the first aid provider and the person on harm by a breathing mask is with a one-way-valve.



Art. no. 5155

## PLUM BURN GEL

Gel dressing that rapidly cools and effectively relief pain for burns. The cooling effect reduces pain and swelling and prevents the burn from going deeper into the skin.



Art. no. 5156

## PLUM QUICKCOOL

Burn gel for first aid in case of minor burns and scalding. Cooling and pain-relieving effect. Contains tea tree oil. Can be used individually and for refilling the first aid stations.



Art. no. 5150

## PLUM QUICKCLEAN

Wound cleansing wipes for cleaning minor wounds and cuts. Each wipe is packed separately and contains sterile water. Can be used individually and for refilling first aid stations.



Art. no. 5151

## PLUM QUICKSTOP

Wound dressing set for minor as well as major bleeding. Sterile dressing with integrated compression. Elastic, non-adhesive and skin-friendly. For general use as pressure, support and protection dressing.



Art. no. 5152

## PLUM WOUND AND EYEWASH SPRAY

A 0.9% phosphate buffered sodium chloride solution, that removes foreign bodies, blood, clot and dead tissue and gives immediate care to minor eye accidents or cuts and scrapes on the skin. 200 ml.



Art. no. 4556

## PLUM QUICKFIX ELASTIC

Elastic and comfortable textile plasters that adjust to the movements of the skin. The plasters are skin-friendly and allow the skin to breathe.



Art. no. 5512  
Art. no. long 5508

## PLUM QUICKFIX DETECTABLE

Blue conspicuous plasters with a built-in metal foil that can be traced by a metal detector. Made from a breathable material.



Art. no. 5513  
Art. no. long 5509

## PLUM QUICKFIX WATER RESISTANT

Water resistant plasters that are suitable for working environments with a high level of humidity. The plasters are made from a protective PE material.



Art. no. 5511

## PLUM QUICKFIX ALOE VERA

Aloe vera contained plasters with soothing, moisturizing and calming effect. Creates the best conditions for natural wound healing and prevents the formation of scabs.



Art. no. 5514

## PLUM QUICKFIX ALU

Air-permeable and comfortable plaster with micronized aluminum pad. Anti-adhering layer in aluminum and PE with high anti-adhering performance.



Art. no. 5515

## PLUM QUICKFIX BLOOD STOPPER

Hemostatic effect with alginate fibres between the nonwoven fibres of the wound pad. Creates artificial clot inside the pad, which blocks bleeding in 1-2 minutes and maintains optimal moist conditions.



Art. no. 5516

## PLUM EYEWASH 0.9% SODIUM CHLORIDE QUICKRINSE

5 eyewash ampoules with 20 ml 0.9% sterile sodium chloride. For use individually or for refilling Plum Pull1Aid®, QuickSafe®, Plum QuickFix&Rinse and Plum QuickRinse dispensers.



Art. no. 5160