Source File

PLUM QUICKFIX IS A UNIQUE DISPENSER:

- Immediate and visible access to plasters and eyewash ampoules
- Each plaster is ready for use with one pull
- Apply with one hand
- Easily refilled as required
- Easy to see when it is empty
- Quickly stops minor bleeding
- Small dimensions easy to fit
- Includes a shelf for storing small, relevant objects
- Can be locked so that it can only be opened with a key

Your quick help for minor everyday accidents

Mount the plaster dispenser on the wall and save time. Put on the plaster - and back to work!

The Plum QuickFix dispenser is developed to fit many different workplaces and can be filled with plasters and ampoules of eyewash in a variety of mixes to meet the needs of every work situation.

The green cap protects the content of the dispenser from dust and dirt. The dispenser has an inside shelf for smaller objects and can be secured with lock and key.

It is designed for mounting on any stationary workplace and can be used in any location.

Dimensions: h: 135 x w: 230 x d: 32 mm. Packaging: 10 pcs.

SCAN & WATCH FREDDY'S SIMPLE GUIDE

QuickFix

plum

7 DTHER



info@plum.eu +45 6916 9600

Choose your Plum QuickFix plaster dispenser

- to meet your specific needs



PLUM QUICKFIX WATER RESISTANT

Dispenser with 90 water resistant plasters suitable for environments with a high level of humidity.

Packaging:	10
Article no:	5501
Refill article no:	5511

PLUM QUICKFIX ELASTIC

the skin.

Packaging:

Refill article no:

Article no:

A

F

Dispenser with 90 elastic plasters

that adjust to the movement of



Quick Fix	
plum	
888	888

ix	PLUM QUICKFIX DETECT Dispenser with 90 blue plaste with built-in metal surface whi be traced by a metal detector	rs ch can
	Packaging:	10

Dispenser with 90 blue plasters
vith built-in metal surface which can
be traced by a metal detector.

10
5503
5513

PLUM QUICKFIX ALOE VERA Dispenser with 90 aloe vera plasters, which soothes and moisturizes.

Packaging:	10
Article no:	5506
Refill article no:	5514

PLUM QUICKFIX ALU

the aluminum layer.

Packaging:

Article no:

Refill article no:

Dispenser with 90 alu plasters

with anti-adhering effect due to





PLUM QUICKFIX
BLOOD STOPPER

Dispenser with 90 blood stopper plasters with hemostatic effect.

Packaging:	10
Article no:	5510
Refill article no:	5516



PLUM QUICKFIX ELASTIC & WATER RESISTANT

Dispenser with 45 elastic plasters and 45 water resistant plasters.

ackaging:	10
rticle no:	5507
lefill article no:	5512 & 5511

PLUM QUICKFIX ELASTIC LONG

Dispenser with 60 long elastic plasters that can be wrapped around the finger several times.

Packaging:	10
Article no:	5528
Refill article no:	5508

PLUM QUICKFIX DETECTABLE LONG

Dispenser with 60 long blue plasters with built-in metal surface which can be traced by a metal detector.

Packaging:	10
Article no:	5529
Refill article no:	5509

PLUM QUICKFIX&RINSE

Dispenser with 45 elastic plasters and 5 ampoules with 20 ml 0.9% sterile sodium chloride.

Packaging:	10
Article no:	4631
Refill article no:	5512 & 5160

PLUM QUICKRINSE

Dispenser with 10 eyewash ampoules with 20 ml 0.9% sterile sodium chloride.

Packaging:	10
Article no:	4632
Refill article no:	5160

PLUM QUICKFIX

EMPTY DISPENSER Dispenser without refills. Fill with any Plum QuickFix Plum Quick-Rinse refills of your choice.

Packaging:	10
Article no:	5500



10

5505

5515







QuickFix&Rinse	

10 5502 5512

P А R

DuickRinse

Quick**Fix**

Plum QuickFix plaster refills



PLUM QUICKFIX ELASTIC

45 elastic and comfortable textile plasters that adjust to the movements of the skin. The plasters are skin-friendly and allow the skin to breathe.



PLUM QUICKFIX DETECTABLE

45 blue conspicuous textile plasters with a built-in metal foil that can be traced by a metal detector. Particularly suitable for the food industry. The plasters are made from a breathable material.



PLUM QUICKFIX WATER RESISTANT 45 water resistant plasters that are particularly suitable for working environments with a high level of humidity. The plasters are made from a protective PE material.

REFILL Display carton: Packaging: Article no: Dimensions:	6 x 45 8 x 6 x 45 5512 I: 72 x w: 25 mm
DISPENSER Article no:	5502



REFILL Display carton: Packaging: Article no: Dimensions:	6 x 45 8 x 6 x 45 5511 I: 72 x w: 25 mm
DISPENSER Article no:	5501



PLUM QUICKFIX LONG ELASTIC AND DETECTABLE

30 textile plasters with extra length that can be wrapped around the finger several times. Choose between elastic plasters or blue detectable plasters.





PLUM QUICKFIX ELASTIC LONG 6 x 30

Display carton:	
Packaging:	
Article no:	
Dimensions:	l: 120 x

DISPENSER Article no:

5528

5508

8 x 6 x 30

w: 20 mm



PLUM DETECTABLE LONG	
Display carton:	6 x 30
Packaging:	8 x 6 x 30
Article no:	5509
Dimensions:	l: 120 x w: 20 mm
DISPENSER	

Article no:

5529



www.plum.eu

Plum QuickFix plaster refills




PLUM QUICKFIX ALOE VERA

The aloe vera contained in the pad has well known soothing, moisturizing and calming properties. It creates the best conditions for the wound environment helping the natural healing process. Prevents the formation of scabs. The material of the plaster is textile.

REFILL	
Display carton:	6 x 45
Packaging:	8 x 6 x 45
Article no:	5514
Dimensions:	l: 72 x w: 25 mm
DISPENSER	
Article no:	5506





Plaster with micronized aluminum pad. The anti-adhering layer in aluminmum and PE has a higher anti-adhering performance compared to a standard PE layer. The plaster is air-permeable and comfortable.

REFILL	
Display carton:	6 x 45
Packaging:	8 x 6 x 45
Article no:	5515
Dimensions:	l: 72 x w: 25 mm
DISPENSER	
Article no:	5505



PLUM QUICKFIX BLOOD STOPPER Plasters with hemostatic effect. The presence of alginate fibres between the nonwoven fibres of the wound pad creates an artificial clot (gel) inside the pad, which blocks the bleeding in 1-2 minutes and maintains optimal moist conditions.

REFILL Display carton: Packaging: Article no: Dimensions:	6 x 45 8 x 6 x 45 5516 I: 72 x w: 25 mm
DISPENSER Article no:	5510





PLUM EYEWASH 0.9% SODIUM CHLO-RIDE QUICKRINSE - 5 X 20 ML 5 eyewash ampoules with 20 ml 0.9% sterile sodium chloride. For use individually or for refilling QuickSafe®, Plum QuickFix-&Rinse and Plum QuickRinse dispensers.

Content:	5 x 20 ml
Packaging:	4 x 5 x 20 ml
Article no:	5160



