



Fast and easy first aid within seconds - for skin and eyes!

Plum Wound and Eyewash Spray is a portable, easy to use, saline-based and pain-free spray. It works at any angle, providing a balanced flow rate of sterile 0.9 % saline solution, strong enough for washing away clot and foreign bodies from wounds, yet soft enough for washing the eyes.

PRODUCT FACTS & BENEFITS

- Sterile until empty - allows repeated use
- Quick and easy single hand operation
- Safe rinsing
- Dust cap to protect the spray nozzle from contamination
- Easy to transport, classed as a non-flammable aerosol
- Up to 100% product emptying
- The can is fully recyclable
- Hypoallergenic and sterile solution without preservatives
- Shelf-life: 3.5 years

AVAILABLE SIZES

Art. no.	45530	4554	4556 - buffer	4123 - wall bracket
Size	50 ml	250 ml	200 ml	for 200 ml + 250 ml
Packaging	24 x 50 ml	12 x 250 ml	12 x 200 ml	24 pcs



Art.no. 4554



Art.no. 455304



Art.no.
4556 + 4123

The product is CE marked according to the Guidelines for Classification of Medical Device Directive (MDD 93/42/EEC).

info@plum.eu
+45 6916 9600