

Plum iBox

- designed for your safety

Delivered with wall bracket and wall mounting set




Quick and efficient overview of eyewash bottles and expiry dates

Release from wall bracket and quickly bring to injured person

2026 08

Quickly accessible eyewash solutions

2026 08



Delivered with an eyewash plan

Strong and durable 05 material (EPP)

Spacious design lets you grab bottles fast and easily

Easily opens - mounted on the wall or standing on the floor



Dust-proof and portable eyewash box

Plum iBox is balancing key features in a modern design and is built in a strong and durable material. See-through window lets you keep track of expiry dates immediately. Release Plum iBox from the wall and quickly bring the eyewash to the injured person.





PRODUCT FACTS & BENEFITS

- Clearly visible eyewash bottles
- Quick and efficient overview of expiry dates and product condition
- Dust-proof box IP55
- Grab from the wall bracket and quickly bring to injured person
- Spacious design lets you grab bottles fast and easy
- Easily opens mounted on the wall or standing on the floor
- Strong and durable material (EPP)
- Delivered with wall bracket, eyewash plan and wall mounting set

AVAILABLE COMBINATIONS		Empty	Plum Eyewash pH Neutral 4.9% Phosphate		Plum Eyewash 0.9% Sodium Chloride		Plum Wound and Eyewash Spray
	Art. no.		200 ml	500 ml DUO	500 ml	1000 ml DUO	250 ml
Plum iBox 1	4900	Χ					
Plum iBox 1	4901					1	
Plum iBox 1	4902				1		
Plum iBox 2	4920	Χ					
Plum iBox 2	4921		1		1		
Plum iBox 2	4922				2		
Plum iBox 2	4923			1		1	
Plum iBox 2	4924					2	
Plum iBox 3	4930	Χ					
Plum iBox 3	4931		1		1	1	1

Install your Plum iBox less than 20 meters from chemical hazard. Release the Plum iBox from the wall and quickly bring it to the injured person. Every second counts!



