# **Technical Data Sheet**

# Plum QuickSoft Bandage (blue, beige, black)

# 1. Identification of the product and the company

#### 1.1. Product identifier and intended use.

Trade Name: Plum QuickSoft Bandage

**Intended Use:** Intended to protect fingers and other joints.

**Expiry date:** Printed on the paper roll.

### 1.2. Details of the supplier of the data sheet

Manufacturer

Company: Plum Safety ApS
Address: Mandelalleen 1

 Zip code:
 5610

 City:
 Assens

 Country
 DENMARK

 E-mail
 info@plum.eu

 Phone
 +45 6916 9600

 Fax
 +45 6916 9630

#### 1.3 Legislative information

General information: As the products are class I medical devices according to Regulation (EU) 2017/745 of

the European Parliament and of the Council of 5 April 2017 on medical devices concerning medical devices. No notified body is involved. No Safety Data Sheet is

required for medical devices.

## 2. General information about the product

#### Physical properties and composition:

	Blue	Black	Beige
Material	60 % Polyester	60 % Polyester	60 % Polyester Polyurethane
	Polyurethane Foam	Polyurethane Foam	Foam
	40% synthetic latex	40% synthetic latex	40% synthetic latex
Healthy Effects	May cause skin irritation	May cause skin irritation	May cause skin irritation.
Physical and Chemical	Flash point 470° C	Flash point 470° C	Flash point 470° C
Properties			
Primary Packing (Pouch)	Packed with 2 layers of	Packed with 2 layers of	Packed with 2 layers of paper
	paper	paper	
Dimensions	5.5 cm x 5 m	5.5 cm x 5 m	5.5 cm x 5 m
Special feature	Latex free	Latex free	Latex free
	Easy tear	Easy tear	Easy tear
	Cohesive/self-adhesive	Cohesive/self-adhesive	Cohesive/self-adhesive

# **Technical Data Sheet**

# Plum QuickSoft Bandage (blue, beige, black)

## 3. Side effects and malfunction

Side effects and malfunction:

Since the product is not sterile it should not be used on large wounds and skin disease without a suitable sterile dressing. This product has a certain elasticity; therefore, it requires appropriate pressure and several layers when binding or fixing. Not too tight or too loose. Too tight will cause discomfort or adverse reactions to the user, such as local ischaemia. Too loose, will not hold the dressing in place or wrap it effectively.

The product is designed to protect fingers and other joints. If side effects or malfunctions are experienced, please report these to Plum Safety ApS, <u>info@plum.eu</u> or phone +45 6471 2112, so that necessary remedial measures can be initiated. Serious adverse effects should always be reported to a relevant and competent authority.

## 4. Storage

Storage:







Keep dry

Keep away from sunlight

Storage temperature

**Shelf life:** 

3 years from production date.

### 5. Other information

#### Version history and indication of changes

Version	Revision date	Responsible	Changes
1.0	14-02-2024	LBH	

Other information Additional information can be obtained by contacting info@plum.eu or by telephone +45 6471 2112.

TDS prepared by:

Company: Plum Safety ApS
Address: Mandelalleen 1

**Zip code:** 5610

City

Assens

 Country:
 DENMARK

 E-mail:
 info@plum.eu

 Phone:
 +45 6916 9600

 Fax:
 +45 6916 9630

Document Language: English