



SAFETY AT WORK

Plum Pull1Aid®



info@plum.eu
+45 6916 9600

plum
Safety



Plum Pull1Aid®

Quick, easy and safe first aid

Easy access to first aid is the key to increase safety at work

Plum Pull1Aid® is a wall mounted first aid station that fits most workplaces. A flexible solution that gives a quick overview of the available first aid - and makes it easy to always have the right products in the right place.

Plum Pull1Aid® allows you to combine first aid and wound care products that suit the needs at your workplace.

PRODUCT FACTS & BENEFITS

- ✓ Clear cover - hygienic and protects the products from dust
- ✓ Practical and functional station
- ✓ High quality products
- ✓ Simple to use with clear instructions
- ✓ Easy and quick to refill
- ✓ Immediate and visible access to plasters, ready to use with a pull
- ✓ Wound care products are quick to pull out of the stations

Plum Pull1Aid® is clear, fast and simple first aid so you can act immediately if an accident happens. Expand your Plum Pull1Aid® station for increased safety for wounds and eyes by mounting a QuickSafe® station next to it.

Plum Pull1Aid® Station



CONTENTS:

1 x Plum QuickCool burn gel	5150
1 x Plum QuickClean wound cleanser	5151
1 x Plum QuickStop wound dressing	5152
2 x Plum Blood Stopper wound dressing	5154
1 x Plum Blood Stopper Mini wound dressing	5153
1 x Plum BurnGel dressing for burns	5156
1 x Plum Protect Kit respiratory aid	5155
1 x Plum QuickFix Elastic refill	5512
1 x Plum Eyewash 0.9% Sodium Chloride QuickRinse	5160

ART. NO.:

Packaging:	1
Article no:	5183
Dimensions:	h: 430 x w: 253 x d: 92 mm

Plum Pull1Aid® Refills

PLUM BLOOD STOOPER MINI

A smaller, sterile compress dressing with four functions. Use as pressure pad for heavy bleeding, protection for superficial wounds, burn dressing with wound-protecting surface or as temporary support.



Art. no. 5153

PLUM BLOOD STOPPER

A sterile compress dressing with four functions. Use as pressure pad for heavy bleeding, protection for superficial wounds, burn dressing with wound-protecting surface or as temporary support.



Art. no. 5154

PLUM PROTECT KIT

Protecting aid for mouth-to-mouth resuscitation situations. Limits contact between the first aid provider and the person on harm by a breathing mask is with a one-way-valve.



Art. no. 5155

PLUM BURN GEL

Gel dressing that rapidly cools and effectively relief pain for burns. The cooling effect reduces pain and swelling and prevents the burn from going deeper into the skin.



Art. no. 5156

PLUM QUICKCOOL

Burn gel for first aid in case of minor burns and scalding. Cooling and pain-relieving effect. Contains tea tree oil. Can be used individually and for refilling the first aid stations.



Art. no. 5150

PLUM QUICKCLEAN

Wound cleansing wipes for cleaning minor wounds and cuts. Each wipe is packed separately and contains sterile water. Can be used individually and for refilling first aid stations.



Art. no. 5151

PLUM QUICKSTOP

Wound dressing set for minor as well as major bleeding. Sterile dressing with integrated compression. Elastic, non-adhesive and skin-friendly. For general use as pressure, support and protection dressing.



Art. no. 5152

PLUM QUICKFIX BLOOD STOPPER

Hemostatic effect with alginate fibres between the nonwoven fibres of the wound pad. Creates artificial clot inside the pad, which blocks bleeding in 1-2 minutes and maintains optimal moist conditions.



Art. no. 5516

PLUM QUICKFIX ELASTIC

Elastic and comfortable textile plasters that adjust to the movements of the skin. The plasters are skin-friendly and allow the skin to breathe.



Art. no. 5512
Art. no. long 5508

PLUM QUICKFIX DETECTABLE

Blue conspicuous plasters with a built-in metal foil that can be traced by a metal detector. Made from a breathable material.



Art. no. 5513
Art. no. long 5509

PLUM QUICKFIX WATER RESISTANT

Water resistant plasters that are suitable for working environments with a high level of humidity. The plasters are made from a protective PE material.



Art. no. 5511

PLUM QUICKFIX ALOE VERA

Aloe vera contained plasters with soothing, moisturizing and calming effect. Creates the best conditions for natural wound healing and prevents the formation of scabs.



Art. no. 5514

PLUM QUICKFIX ALU

Air-permeable and comfortable plaster with micronized aluminum pad. Anti-adhering layer in aluminum and PE with high anti-adhering performance.



Art. no. 5515

PLUM EYEWASH 0.9% SODIUM CHLORIDE QUICKRINSE

5 eyewash ampoules with 20 ml 0.9% sterile sodium chloride. For use individually or for refilling Plum Pull1Aid®, QuickSafe®, Plum QuickFix&Rinse and Plum QuickRinse dispensers.



Art. no. 5160

PLUM PULL1AID - EMPTY STATION

Practical and functional station with clear cover - hygienic and protects the products from dust.



Art. no. 5182