

Plum Wound and Eyewash 0.9% Sodium Chloride QuickRinse



SAFETY AT WORK

Plum QuickRinse - for smaller accidents.

Plum Wound and Eyewash 0.9% Sodium Chloride QuickRinse is a sterile first aid wound and eye irrigation.

The product can be used for washing away:

- Dust
- Dirt
- Foreign bodies

Plum Wound and Eyewash 0.9% Sodium Chloride QuickRinse ampoules are designed specifically for rinsing away small foreign particles that can be removed with just a small amount of fluid. A key advantage of these ampoules is their convenience – you only need to open a small, single-use container rather than a larger bottle.



PRODUCT FACTS & BENEFITS

- Effectively removes dust, dirt, and contamination
- Easy to handle - easily applied by the injured person or with assistance
- Shelf life of 3.5 years for sealed ampoules



This is a sterile 0.9 % sodium chloride solution that eases discomfort and prevents subsequent complications caused by foreign bodies.

A gentle solution


The 0.9% sodium chloride solution matches the eye's natural saline levels, which makes it gentle and comfortable for the user, reducing any unpleasant sensation during rinsing.

INSTRUCTIONS FOR USE

Plum Wound and Eyewash 0.9% Sodium Chloride QuickRinse is easy to use. Follow these simple steps:



Step 1: Twist the ampoule.



Step 2: Rinse eye or wound.



Step 3: Throw away the empty ampoule.



PLUM WOUND AND EYEWASH 0.9% SODIUM CHLORIDE QUICKRINSE - 5 X 20 ML

5 ampoules with 20 ml 0.9% sterile sodium chloride. For use individually or for refilling Plum QuickSafe®, Plum QuickFix&Rinse and Plum QuickRinse dispensers.

Content: 5 x 20 ml
Packaging: 4 x 5 x 20 ml
Article no: 5160



PLUM QUICKRINSE

Dispenser with 10 ampoules with 20 ml 0.9% sterile sodium chloride.

Packaging: 10
Article no: 4632
Refill article no: 5160



PLUM QUICKFIX&RINSE

Dispenser with 45 elastic plasters and 5 ampoules with 20 ml 0.9% sterile sodium chloride.

Packaging: 10
Article no: 4631
Refill article no: 5512 & 5160



PLUM QUICKFIX EMPTY DISPENSER

Dispenser without refills. Fill with any Plum QuickFix or Plum QuickRinse refills of your choice.

Packaging: 10
Article no: 5500

info@plum.eu
 +45 6916 9600

