



SAFETY AT WORK

Plum iBox

– designed for your safety


Delivered with wall bracket and wall mounting set



Quick and efficient overview of wound and eyewash bottles and expiry dates

Release from wall bracket and quickly bring to injured person

Quickly accessible wound and eyewash solutions



Delivered with a wound and eyewash plan

Strong and durable material (EPP)




Spacious design lets you grab bottles fast and easily

Easily opens - mounted on the wall or standing on the floor

Dust-proof and portable box

Plum iBox is balancing key features in a modern design and is built in a strong and durable material. See-through window lets you keep track of expiry dates immediately. Release Plum iBox from the wall and quickly bring the eyewash to the injured person.



PRODUCT FACTS & BENEFITS

- Clearly visible bottles
- Quick and efficient overview of expiry dates and product condition
- Dust-proof box - IP55
- Grab from the wall bracket and quickly bring to injured person
- Spacious design lets you grab bottles fast and easy
- Easily opens - mounted on the wall or standing on the floor
- Strong and durable material (EPP)
- Delivered with wall bracket, wound and eyewash plan and wall mounting set

AVAILABLE COMBINATIONS

	Art. no.	Empty	Plum Wound and Eyewash pH Neutral 4.9% Phosphate		Plum Wound and Eyewash 0.9% Sodium Chloride		Plum Wound and Eyewash Spray
			200 ml	500 ml DUO	500 ml	1000 ml DUO	250 ml
Plum iBox 1	4900	X					
Plum iBox 1	4901					1	
Plum iBox 1	4902				1		
Plum iBox 2	4920	X					
Plum iBox 2	4921		1		1		
Plum iBox 2	4922				2		
Plum iBox 2	4923			1		1	
Plum iBox 2	4924					2	
Plum iBox 3	4930	X					
Plum iBox 3	4931		1		1	1	1

Install your Plum iBox less than 20 meters from chemical hazard. Release the Plum iBox from the wall and quickly bring it to the injured person. Every second counts!

info@plum.eu
+45 6916 9600