

Technical Data Sheet

Plum QuickSoft Cohesive Bandage (blue, beige, black)

1. Identification of the product and the company

1.1. Product identifier and intended use.

Trade Name:	Plum QuickSoft Cohesive Bandage
Intended Use:	Intended to protect fingers and other joints.
Expiry date:	Printed on the paper roll.

1.2. Details of the supplier of the data sheet

Manufacturer	
Company:	Plum Safety ApS
Address:	Mandelalleen 1
Zip code:	5610
City:	Assens
Country	DENMARK
E-mail	info@plum.eu
Phone	+45 6916 9600
Fax	+45 6916 9630

1.3 Legislative information

General information:	As the products are class I medical devices according to <i>Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices</i> concerning medical devices. No notified body is involved. No Safety Data Sheet is required for medical devices.
-----------------------------	---

2. General information about the product

Physical properties and composition:

	Blue	Black	Beige
Material	60 % Polyester Polyurethane Foam 40% synthetic latex	60 % Polyester Polyurethane Foam 40% synthetic latex	60 % Polyester Polyurethane Foam 40% synthetic latex
Healthy Effects	May cause skin irritation	May cause skin irritation	May cause skin irritation.
Physical and Chemical Properties	Flash point 470 ⁰ C	Flash point 470 ⁰ C	Flash point 470 ⁰ C
Primary Packing (Pouch)	Packed with 2 layers of paper	Packed with 2 layers of paper	Packed with 2 layers of paper
Dimensions	5.5 cm x 5 m	5.5 cm x 5 m	5.5 cm x 5 m
Special feature	Latex free Easy tear Cohesive/self-adhesive	Latex free Easy tear Cohesive/self-adhesive	Latex free Easy tear Cohesive/self-adhesive

Technical Data Sheet

Plum QuickSoft Cohesive Bandage (blue, beige, black)

3. Side effects and malfunction

Side effects and malfunction:

Since the product is not sterile it should not be used on large wounds and skin disease without a suitable sterile dressing. This product has a certain elasticity; therefore, it requires appropriate pressure and several layers when binding or fixing. Not too tight or too loose. Too tight will cause discomfort or adverse reactions to the user, such as local ischaemia. Too loose, will not hold the dressing in place or wrap it effectively.

The product is designed to protect fingers and other joints. If side effects or malfunctions are experienced, please report these to Plum Safety ApS, info@plum.eu or phone +45 6471 2112, so that necessary remedial measures can be initiated. Serious adverse effects should always be reported to a relevant and competent authority.

4. Storage


Storage:



Storage temperature



Keep dry



Keep away from sunlight

Shelf life:

3 years from production date.

5. Other information

Version history and indication of changes

Version	Revision date	Responsible	Changes
1.0	14-02-2024	LBH	
1.1	18-12-2024	MAW	1.1 + header

Other information

Additional information can be obtained by contacting info@plum.eu or by telephone +45 6471 2112.

TDS prepared by:

Company: Plum Safety ApS
Address: Mandelalleen 1
Zip code: 5610
City:

Assens

Country: DENMARK
E-mail: info@plum.eu
Phone: +45 6916 9600
Fax: +45 6916 9630

Document Language: English