

Technical Data Sheet

Pull1Aid® Protect Kit

1. Identification of the product and the company

1.1. Product identifier and intended use

Trade Name:	Pull1Aid® Protect Kit
Intended Use:	Plum Pull1Aid® Protect Kit to limit contact between the first aid provider and the injured person.
Expiry date:	Printed on the sealed package, when needed.

1.2. Details of the supplier of the data sheet







Distributor	
Company:	Plum Safety ApS
Address:	Mandelalleen 1
Zip code:	5610
City:	Assens
Country:	Denmark
E-mail:	info@plum.eu
Phone:	+45 6916 9600
Fax:	+45 6916 9630

1.3 Legislative information

General information:	REF 97106K, Respiratory sheet: Are class 1 medical devices according to <i>Council Directive 93/42/EEC of 14 June 1993</i> concerning medical devices. <i>No Safety Data Sheet is required for medical devices.</i> Legal manufacturer is FIRSTAR HEALTHCARE Co. Ltd. (Guangzhou), Rm 901, Building No. 2, Headquarters Center, Tian'an High-tech Ecological Park, Panyu, 511400, Guangzhou, P.R. China Distributor: AKLA AB, Parkgatan 15, 696 33 Askersund, Sweden Commercial contact: Plum Safety ApS, Mandelalleen 1, 5610 Assens, Denmark REF 94416ST, Disinfection wipe: Reference to MSDS REF 2080809001, Protective Gloves: Not classified
----------------------	--

2. General information about the product

Physical properties and composition:

Pull1Aid®Protect Kit. Art no. 5155	REF. 97106K, Respiratory sheet	REF. 94416ST, Disinfection wipe	REF. 2080809001, Protective Gloves (1 pair)
Material	Valve= PP, Sheet=PE	Nonwoven tissue, Cellulose, Viscose, Co-polyester/Polyefine saturated with 70% alcohol solution (alcohol, aqua, glycerine)	Powder-free vinyl gloves, unsterile, transparent, smooth
Primary packing			
Secondary packing			

Technical Data Sheet

Pull1Aid® Protect Kit

Dimensions	-	17 x 22 cm	Size 8-9 (Ø105xL230mm), Thickness 0,05 mm.
Special features	For mouth-to-mouth ventilation	For disinfection in case of first aid.	For hand protection in case of first aid.

3. Side effects and malfunction

Side effects and malfunction:

The products are designed for skin contact and first aid use involving burns. If side effects or malfunctions are experienced, please report these to Plum Safety ApS, info@plum.eu or phone +45 6916 9600, so that necessary remedial measures can be initiated.

Reporting of accidents and incidents. Vigilance. Serious incidents involving products made available on the EU market must be reported immediately to the National Competent Authority. All serious incidents must be reported no later than 15 days after discovery, in the event of a serious threat to public health within 2 days. In the event of death or an unexpected and serious deterioration of a person's state of health, reporting must take place within 10 days.

4. Storage

Shelf life: 5 years from production date

5. Other information

Do not reuse



Version history and indication of changes

Version	Revision date	Responsible	Changes
1.1	08-07-2025	MAW	1.3 + 2
1.0	20-06-2022	CL	

Other information Additional information can be obtained by contacting info@plum.eu or by telephone +45 6916 9600.

TDS prepared by:

Company: Plum Safety ApS
Address: Mandelalleen 1
Zip code: 5610
City: Assens
Country: DENMARK
E-mail: info@plum.eu
Phone: +45 6916 9600
Fax: +45 6916 9630

Document Language: English