

Technical Data Sheet

Plum QuickStop

1. Identification of the product and the company

1.1 Product identifier and intended use

Trade Name:	Plum QuickStop.
Intended Use:	Plum QuickStop is a sterile first aid product that is intended to be used as a mechanical barrier and compression.
Expiry date:	Printed on the sealed package.

1.2 Details of the supplier of the data sheet

Distributor	
Company:	Plum Safety ApS
Address:	Mandelalleen 1
Zip code:	5610
City:	Assens
Country:	Denmark
E-mail:	info@plum.eu
Phone:	+45 6916 9600
Fax:	+45 6916 9630

1.3 Legislative information

General information: As the products are class Is sterile medical devices according to *Council Directive 93/42/EEC of 14 June 1993 and MDR2017/745 art. 120* concerning medical devices, they are certified by the Notified Body, Intertek Medical Notified Body AB, Kista, Sweden.

No Safety Data Sheet is required for medical devices.

Plum Safety ApS is certified according to *ISO 13485:2016 Medical devices – Quality management systems – Requirements for regulatory purposes*.

AKLA AB, Parkgatan 15, 696 33 Askersund, Sweden is the legal manufacturer of Plum QuickStop.

2. General information about the product


Material requirements: The materials in these products have been judged biocompatible according to the requirements of *ISO 10993-1:2020 Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process*.

Physical properties and composition

	Plum QuickStop Art. No. 5152 (REF 92111/92109)
Material	<i>Nonwoven dressing, large, in CELLOSOFT-material: 70% Viscose 30% Polyethylene Non-adhesive, stitched on elastic bandage 4m x 6cm Extra bandage: 3m x 6cm Bandage: Polyamid 60%/Viscose 40%</i> <i>Nonwoven dressing, small, in CELLOSOFT-material, 3-layer combination of: 100% LDPE-HDPE 30g/m² micro-perforated. 100% Viscose 200 g/m² +/- 10% 100% LDPE-HDPE 30g/m² micro-perforated Non-adhesive, stitched on elastic bandage 2m x 6cm</i>
Primary Packing (Pouch)	<i>Paper Plastic: Polyamide 50/Polyethylene COEX 68</i>

Technical Data Sheet

Plum QuickStop

Secondary Packing (Box)	
Absorption	≥ 6g/g
Dimensions	1 pcs 17x17 cm 2 pcs 7.5x12 cm
Special features	Suitable as a pressure bandage and protective dressing for hand, finger, toes

3. Side effects and malfunction

Side effects and malfunction:

The products are designed for skin contact and first aid use involving wounds. If side effects or malfunctions are experienced, please report these to Plum Safety ApS, info@plum.eu or phone +45 6916 9600, so that necessary remedial measures can be initiated.

Reporting of accidents and incidents. Vigilance. Serious incidents involving products made available on the EU market must be reported immediately to the National Competent Authority. All serious incidents must be reported no later than 15 days after discovery, in the event of a serious threat to public health within 2 days. In the event of death or an unexpected and serious deterioration of a person's state of health, reporting must take place within 10 days.

4. Storage

Shelf life: 5 years from production date

5. Other information

Do not reuse



Do not use if package is damaged



Sterilized using irradiation



Version history and indication of changes

Version	Revision date	Responsible	Changes
1.0	14-04-2026	MAW	-

Other information Additional information can be obtained by contacting info@plum.eu or by telephone +45 6916 9600.

TDS prepared by:

Company: Plum Safety ApS
Address: Mandelalleen 1
Zip code: 5610
City: Assens
Country: Denmark
E-mail: info@plum.eu
Phone: +45 6916 9600
Fax: +45 6916 9630

Document Language: English